
Cryoablation Console Operator's Manual

**Model:
Cryoablation Console**

Confidentiality Statement

The information contained in this document is the property of Adagio Medical, Inc. It is therefore provided to you in confidence. It is understood that this information will not be disclosed to others without written authorization from the Adagio Medical Inc., Clinical and/or Regulatory Affairs.

WARRANTY.....	3
CLASSIFICATION	3
TERMS AND DEFINITION	3
1. INTRODUCTION	4
1.1. Intended Use/ Indication For Use	4
1.2. Intended Users of the Device	4
1.3. Intended Patient Population.....	4
1.4. Clinical Benefits	4
1.5. Device Description	4
1.5.1. Overview	4
1.5.2. Console Operation.....	4
1.5.3. GUI.....	5
1.5.4. Console Connections.....	5
1.6. Important Safety Considerations.....	5
1.6.1. Warnings and Precautions	5
1.6.2. Contraindications.....	6
1.6.3. Limitations	6
1.6.4. Adverse Events.....	6
1.7. Training.....	6
1.8. Device Compatibility.....	6
2. PRE-OPERATING INSTRUCTIONS	9
2.1. Unpacking and Setup.....	9
2.1.1. Console	9
2.1.2. Hose Assemblies.....	9
2.1.3. Nitrogen Tanks	9
2.1.4. Handling / Transporting Compressed Gas Cylinder	9
2.1.5. Storage of Compressed Gas Cylinders.....	9
2.1.6. Use of Personal Protective Equipment (PPE).....	10
2.1.7. Cryogen Tank.....	10
2.1.8. Handling Liquid Nitrogen	10
2.1.9. Note.....	10
2.2. Transporting, Positioning and Locking the Console Wheels	11
3. OPERATING INSTRUCTIONS - ADAGIO MEDICAL CRYOABLATION CONSOLE (AF).....	12
3.1. Procedure Overview	12
3.2. Pre-Cryoablation Steps	12
3.2.1. "ABOUT TAB"	16
3.3. "PREPARE CONSOLE" FOR A PROCEDURE.....	16
3.4. "NEW PATIENT" AND CONNECT CATHETER.....	23
3.5. Cryoablation Steps	33
3.5.1. Procedure Screen Numerical Values.....	33
3.5.2. Initiating a Freeze.....	33
3.5.3. Terminating Freeze	37
3.6. Post-Cryoablation Steps.....	38
3.6.1. Disconnecting a Catheter.....	38
3.6.2. Console Shutdown	40
4. OPERATING INSTRUCTIONS - ADAGIO MEDICAL CRYOABLATION CONSOLE (VT)	42
4.1. Procedure Overview	42
4.2. Pre-Cryoablation Steps	42
4.2.1. "ABOUT TAB"	45

4.3. "PREPARE CONSOLE" for a Procedure	46
4.4. "NEW PATIENT" and Connect Catheter	53
4.5. Cryoablation Steps	60
4.5.1. Procedure Screen Numerical Values	60
4.5.2. Initiating a Freeze	60
4.5.3. Terminating Freeze	63
4.6. Post-Cryoablation Steps	64
4.6.1. Disconnecting a Catheter	64
4.6.2. Console Shutdown	66
5. SPECIFICATIONS.....	68
5.1. General Mechanical Specifications.....	68
5.2. Electrical Specifications	68
5.3. Operational Specifications and Environment.....	68
5.4. Transport Environment	68
5.5. Technical Data	68
5.6. Console Configuration Kit	69
6. MANUFACTURER'S DECLARATION	70
6.1. Manufacturer's declaration – Electromagnetic Emission (Table 2, IEC 60601-1-2:2014, 4.1 Edition)	70
6.2. Manufacturer's declaration – Electromagnetic Immunity I (Table 4, 5, and 7 IEC 60601-1-2:2014, 4.1 Edition)	70
6.3. Manufacturer's declaration – Electromagnetic Immunity II (Table 5, 7, and 8, IEC 60601-1-2:2014, 4.1 Edition)	71
6.4. Manufacturer's declaration – Recommended Separation Distances Between Portable and Mobile RF-Communications Equipment and the Product (Table 9, IEC 60601-1-2:2014, 4.1 Edition)	71
7. MAINTENANCE AND SERVICING.....	72
7.1. Preventive Maintenance.....	72
7.1.1. Annual Maintenance.....	72
7.1.2. 5-year Service Maintenance	72
7.2. Software Updates	72
7.3. Refrigerant Shipping	72
7.4. Safe Disposal.....	72
8. INFORMATIONAL MESSAGES, ERRORS, AND TROUBLESHOOTING	73
8.1. Informational Messages	73
8.2. Error Pop-ups.....	74
8.3. Safe Abort Screens	75
8.4. Working with System Informational Messages and Error Pop-ups	76
8.5. Troubleshooting the Monitor	76
8.6. Console Failure during Procedure	76
8.7. Sensor Failure	76
8.8. Catheter Failure	76
8.9. Software Failure.....	76
9. CLEANING.....	77
10. SYMBOLS	77
11. LABELING	79
12. INCIDENT REPORTING	81

WARRANTY

Unless otherwise specified Adagio Medical, Inc. warrants all products of its manufacture at the time of shipment to be free from defects of material and workmanship for a period of twelve (12) months from the date of shipment when owned by the original purchaser. Any product which is believed to be defective, if returned within twelve (12) months after date of shipment by the Company with freight prepaid and found by the Company's inspection to be defective within the terms of this warranty, will be repaired or replaced free of charge and shipped, freight prepaid, to any point in the United States. If inspection by the Company of any such product does not disclose any defect within the terms of this warranty, the Company's regular charges for repairs or replacement and freight shall apply. All consumable and disposable products are guaranteed to be free from defects upon shipment only. The warranty period for batteries is limited to Ninety (90) days from date of shipment.

This warranty specifically excludes and replaces all other express and implied warranties. The Company shall have no liability under any warranty in any amount greater than the Company receives for the sale of the product involved. The use and manner of use of the Company's products shall be the responsibility of the purchaser and the purchaser covenants and agrees to indemnify and save harmless the Company in respect to any loss and damage that may arise through the use by the purchaser, or others, of any of the Company's products.

This warranty is rendered void and Adagio Medical cannot be held liable for conditions resultant therefrom if:

1. Damage to the Console is incurred as a result of modification or mishandling.
2. The customer fails to maintain the Console in a proper manner.
3. The customer uses any parts, accessories or fittings not specified or sold by Adagio Medical, Inc.
4. Sale or Service is performed by a non-certified service/dealer agency or any other unauthorized agency.

All shipping claims must be made within thirty (30) days from date of shipment from Adagio Medical, Inc. otherwise factory will not be liable for claims of missing items. Any item ordered in error and returned to the factory for credit, will be subject to a minimum restocking charge of 15%. Requests for returning items must be made within thirty (30) days of factory shipment date.

Adagio Medical will accept no returned goods without a Returned Goods Authorization number (RGA) obtained from Customer Service Department.

CLASSIFICATION

According to the standard IEC 60601-1 of the International Electrotechnical Commission, *Medical electrical equipment, Part 1: General requirements for safety*, the Adagio Cryoablation Console is classified as follows:

- Type CF, according to the degree of protection against electric shock.
- Equipment is not suitable for use in the presence of a flammable anesthetic mixture with air or with oxygen or nitrous oxide.

TERMS AND DEFINITION

Term	Meaning
Dewar	The Console has a Thermal Reservoir inside, which is used to hold the Liquid Nitrogen and Gas Nitrogen used during procedures. Dewar is the technical name for this Reservoir, which works much like a thermos bottle.
Refrigerant	Liquid Nitrogen (LN2) is an extremely cold fluid (-196°C) that is stored in the console reservoir (Dewar) and is used to chill the Nitrogen gas to extremely low temperatures.
Cryo	From the Greek word for "Cold" meaning extremely low temperatures (-196°C) which the Cryoablation Catheter provides during procedures.
GUI	Acronym for "Graphical User Interface," which allows the user to control and monitor the System.
N2	N2 is the Chemical Formula for Nitrogen. Nitrogen is an element which is found in our atmosphere. It is safe, stable, and noncombustible. The Console uses Liquid Nitrogen and Nitrogen gas in the Cryoablation Catheter during procedures.
NCN (also CN2)	Near-Critical Nitrogen is the pressurized Nitrogen gas which has been chilled to an extremely low temperature (-190°C) through conduction heat transfer from the LN2. In the console closed loop system, the NCN is circulated from the console to the catheter to ablate arrhythmogenic tissues during procedures then returned to the console to be safely bled out to ambient environment.

Table 1: Terms and Definitions

1. INTRODUCTION

1.1. Intended Use/ Indication For Use

Intended Use - The Adagio Medical Inc. Cryoablation Console is intended to be used with the Adagio ablation Catheter(s) for the delivery of cryogenic energy to create transmural, linear, and focal lesions for the necrosis (ablation) of cardiac tissues.

Indication for Use (AF): The Adagio Medical Inc. iCLAS™ Cryoablation System (Catheter, Stylets, and Console) is indicated for the treatment of drug refractory, recurrent, symptomatic, Paroxysmal Atrial Fibrillation (PAF), Persistent Atrial Fibrillation (PsAF), and Atrial Flutter (AFL) by ablation of arrhythmogenic tissue that drives and maintains these arrhythmias.

Indication for Use (VT): The Adagio Medical Inc. VT Cryoablation System (Catheter and Console) is indicated for treatment of monomorphic ventricular tachycardia by ablation of arrhythmogenic tissue that drive and maintain these arrhythmias.

1.2. Intended Users of the Device

The Cryoablation Console is designed for use by Adagio Medical trained field clinical personnel acting under a physician's order.

1.3. Intended Patient Population

The Cryoablation Console is used with the Catheters as a System.

- Intended patient (AF): The Adagio Medical Inc. iCLAS™ Cryoablation System (Catheter, Stylets, and Console) are to be used for treatment of patients suffering from drug refractory, recurrent, symptomatic, paroxysmal atrial fibrillation (PAF), persistent atrial fibrillation (PsAF), or atrial flutter (AFL). Refer to the Contraindications section below for more details.
- Intended patient (VT): The Adagio Medical Inc. VT Cryoablation System (Catheter and Console) is used in patients suffering from recurrent monomorphic ventricular tachycardia. Reference section Contraindication for more details.

1.4. Clinical Benefits

The intended clinical benefit of the Adagio Medical iCLAS Cryoablation System including the Catheter, Stylets, and Console is to prevent the recurrence of atrial tachycardia (AF/AT/AFL episodes). In the CryoCure2 Clinical Study (NCT 02839304), the Catheter and Stylets have demonstrated ability to ablate at various locations in the left and right atria supported by acute success rate of 97% Pulmonary Vein (PV) isolation, 100% posterior wall isolation, 92.3% BDB at cavo-tricuspid isthmus (CTI), and 70% mitral isthmus isolation. At 1-year follow up, 82.6% of the study population including PAF and PsAF patients remained free from AF recurrence after a single cryoablation procedure. Sub-set analysis of the PsAF cohort showed 85.9% of patients did not have AF recurrence at the end of the follow up period.

The intended clinical benefit of the Adagio Medical VT Cryoablation System is to prevent the reoccurrence of sustained VT episodes or appropriate ICD Intervention. In the CryoCureVT Clinical Study (NCT 04893317), acute success was reported at 94% non-inducible VT at the end of the ablation procedure. At 6-month follow up, 63% of the study population remained free from VT recurrence or appropriate ICD intervention after a single cryoablation procedure.

1.5. Device Description

1.5.1. Overview

The Adagio Cryoablation Console is a re-usable, non-sterile, electromechanical device providing the cryoablation energy and control mechanisms. The Console contains an internal reservoir for the Liquid Nitrogen. The Nitrogen is pressurized and delivered to the Cryoablation Catheter creating the cryogenic temperatures (-196°C) to ablate cardiac tissue. The console may be connected to the hospital's nitrogen gas distribution system or to a stand-alone nitrogen gas tank.

The Console is a microcontroller-based device designed to control all operational aspects of the Adagio Cryoablation Catheter such as temperature, pressure, vacuum, nitrogen gas delivery, cycle duration, and the heat exchange. The software will control the delivery of nitrogen gas through control of system valves to achieve target temperatures and pressure. The software will continuously read and monitor system and catheter temperatures and pressures. The software will display ablation parameters and procedure information on a touch screen monitor that is controlled by a PC. To ease transport, the monitor is capable of being folded down, and wheels have been fitted to the console. The wheels can be locked at any time to prevent unintended movement.

The Console is intended to be used with the Adagio ablation Catheter(s) as a system for the delivery of cryogenic energy to create a transmural, linear and focal lesions for the necrosis (ablation) of cardiac tissues, providing greater assurance that the arrhythmogenic substrates are eliminated or isolated, thus restoring a normal cardiac rhythm. The System achieves its intended purpose by creating and controlling the conditioning of Nitrogen to its critical point (CN2) within the Console and subsequently creating an uninterrupted flow of this cooling agent (cryogen) through the Catheter(s).

1.5.2. Console Operation

The Cryoablation Console contains an internal reservoir (Dewar) of Liquid Nitrogen (LN2) refrigerant. This LN2 is used to chill Nitrogen gas to extremely low (cryo) temperature. When Nitrogen gas is chilled to cryo-temperature and pressurized to approximately 500 psi, it becomes Near Critical (NCN) and can circulate through the insulated length of the Catheter and flow into the Cryoablation Element (most distal end). The Cryoablation Element portion is cooled to approximately -190°C and can ablate cardiac tissue. The Catheter is a closed system: after the NCN flows back, the flow is directed into the Console where it is heated, expanded, and safely vented to the atmosphere.

1.5.3. GUI

The Console contains an internal Computer which is used to properly sequence the series of solenoid-control valves that allow the NCN to flow through the Catheter. The Computer functions via a custom software application which interfaces a GUI Touch Screen that allows the Operator to control and monitor the System. All Console operations such as Catheter Connection/Disconnection, Freeze Start/Stop, Settings, and supply levels are sequenced and monitored by the Computer via “push buttons” icons on the GUI.

1.5.4. Console Connections

The Console is designed to operate on external AC power (120, 220, 240 volts 50-60 Hz, factory-set). The console must have an external Nitrogen gas supply of > 800 psi to operate internal pneumatics and maintain sufficient Cryoablation circulation flow through the catheter. Refer to System Specifications and to Console Diagram for more details about connection ports.

1.6. Important Safety Considerations

The Adagio Cryoablation Console has been tested and found to comply with limits for electromagnetic interference and susceptibility as defined by IEC 60601-1-2:2014. If the Cryoablation Console is subjected to Electromagnetic Disturbances, the delivery of cryogenic energy may be impacted and result in the patient not receiving the intended treatment. Refer to Manufacturer's Declaration (Section 6).

The Catheter used on this Console must be disposed of properly.

WARNINGS, CAUTIONS, and NOTES used in this manual have the following significance:

WARNING: Maintenance or operating procedures, techniques, etc., which may result in personal injury or loss of life if not carefully followed.

CAUTION: Maintenance or operating procedures, techniques, etc., which may result in patient harm or damage to equipment if not carefully followed.

NOTE: A maintenance or operating procedure, technique, etc., which is considered essential to emphasize.

The principal WARNING and CAUTION notices to be observed in use of this Console are brought together here for emphasis.

1.6.1. Warnings and Precautions

- **FOR USE BY TRAINED PERSONNEL ONLY** – The equipment contained in this Console is designed for use by trained medical personnel acting under a physician's order.
- **OBSERVE BEST PRACTICE** – The instructions in this manual in no way supersede established medical procedures or staff preference concerning patient care. “Best Practice” as determined by the medical community is to be observed.
- No modification of this equipment is allowed.
- **PRESSURIZED NITROGEN GAS USED IN CONSOLE** – Incorrect operation may result in damage to the device, ineffective treatment, or harm to the patient and user. Do not attempt to disconnect pressurized Nitrogen Gas line while console is in use.
- **SECURE ALL TANKS** – Nitrogen gas tanks must be properly secured. Nitrogen gas is under high pressure.
- **ENSURE CONSOLE IS PROPERLY GROUNDED** – To ensure grounding reliability, connect the power cord to only a properly grounded, 3-wire hospital grade or hospital use outlet of the proper voltage and frequency. DO NOT USE EXTENSION CORDS. Do not connect to console if the power cord shows signs of wear, aging, cracking, or fraying. The use of power cord other than the Adagio provided cord or hospital grade power cable could result in electromagnetic emissions or decrease electromagnetic immunity of the Cryoablation Console and result in improper operation.
- Under any type of sensor failure, the operator must determine the trouble and be extremely attentive to the condition of the patient, if the Console must remain in operation. Never remove a catheter from a patient if ice remains on the Cryoablation Element of the catheter.
- To avoid the risk of electric shock, this equipment must only be connected to a supply main with protective earth.
- **Electromagnetic Disturbances:** Portable and mobile RF communications equipment should be used no closer than 30 cm (12 in) to any part of the Console, including cables, and the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.
- The Cryoablation Console is IPX0 rated and does not provide protection against intrusion of liquids. Do not subject to liquids that could result in damage to the Console or Console electronics and could result in harm to patient or user.
- **SERVICE ONLY BY ADAGIO MEDICAL QUALIFIED PERSONNEL** – The Console should be serviced only by Adagio Medical qualified personnel. Dangerous voltages may be inside and there are no user-serviceable parts.
- Use Console only with Adagio Medical supplied catheters. The catheter is a sterile device and must only be opened and used in the sterile field. After exposure to biological substances, it must be properly disposed.
- The Adagio Esophageal Warming Balloon Catheter is a standalone medical device recommended for use with the Adagio iCLAS™ Cryoablation System to mitigate the risk of a thermal esophageal injury. Do not attempt to operate the Esophageal Warming Balloon Catheter before thoroughly reading the Instructions for Use (Adagio P/N: 108-0053-XXX). As a general safety precaution, the Console should be turned off and AC power disconnected before cleaning.
- Entry of fluids into the Exhaust Fan motor will significantly reduce service life. Clean this area by wiping with a cloth dampened with an appropriate cleaning solution and wipe dry.

- **CAUTION** – Use only Liquid Nitrogen and 99.998% pure Nitrogen Gas (Ultra High Purity) with the Console.
- **CAUTION** – Casters are subject to wear and tear and should be inspected for damage before each transport. Inspect for worn or damaged rubber, bent or loose shafts and tightness of retaining nuts if equipped. Report damaged casters immediately.
- **CAUTION** – Always check the Console is installed in a safe and stable position. Do not position the Console in any orientation other than with the wheels locked and sitting on a stable, flat, horizontal surface.
- **CAUTION** – When moving Cryoablation Console, beware of the Console tipping/over-balancing on an incline of 10 degrees or greater.
- **CAUTION** – Entry of fluids into the Catheter Connection Port can potentially damage the port and cause loss of pressure or vacuum. Never use any type of cleaning fluid on the Catheter Connection Port.
- **CAUTION** – When used with patients who have an implantable medical electronic device (e.g. pacemaker, defibrillator, etc.), be certain the device is not active.

1.6.2. Contraindications

Contraindications (AF): The Adagio Medical Inc. iCLAS™ Cryoablation System (Catheter, Stylets, and Console) is contraindicated as follows:

- In patients with a systemic infection
- In patients with intracardiac thrombus
- In patients with cryoglobulinemia
- In pediatric patients (below the age of 18)
- In pregnant or breastfeeding patients
- In patients with hypercoagulopathy or an inability to tolerate anticoagulation therapy during an electrophysiology procedure
- In the ventricle of the heart where the device may become entrapped in the valve or cordae structures
- In patients with a recent ventriculotomy or atriotomy because this may increase the risk of cardiac perforation or embolic event
- In patients with pulmonary vein stents as the Catheter may dislodge or damage the stent
- In conditions where insertion into or manipulation in the atria is unsafe as this may increase the risk of perforation or systemic embolic event
- In patients with an interatrial baffle or patch as the transseptal puncture could fail to close
- In patients with a contraindication to an invasive electrophysiology procedure where insertion or manipulation of a catheter in the cardiac chambers is deemed unsafe.

Contraindications (VT): The Adagio Medical Inc. VT Cryoablation System (Catheter and Console) is contraindicated as follows:

- In patients with a systemic infection
- In patients with intracardiac thrombus
- In patients with cryoglobulinemia
- In pediatric patients (below the age of 18)
- In pregnant or breastfeeding patients
- In patients with hypercoagulopathy or an inability to tolerate anticoagulation therapy during an electrophysiology procedure
- In patients with an interatrial baffle or patch as the transseptal puncture could fail to close
- In patients with a contraindication to an invasive electrophysiology procedure where insertion or manipulation of a catheter in the cardiac chambers is deemed unsafe.

1.6.3. Limitations

There are no known limitations of the cryoablation console beyond those described in this section (1) including intended use/indication, contraindications and warnings, device compatibility, etc.

1.6.4. Adverse Events

Reference the Instructions for Use for the catheter you are using.

1.7. Training

Prior to use the Operator must read and understand these instructions.

1.8. Device Compatibility

The Cryoablation Console is intended for use with the following compatible Adagio Medical, Inc. catheters only:

- iCLAS (AF4) Cryoablation Catheter
- Esophageal Warming Balloon Catheter
- iCLAS CoreTx (AF3) Cryoablation Catheter
- vCLAS (VT) Cryoablation Catheter

For information about the cryoablation catheter, reference the Instructions for Use.

For information about setting up the Cryoablation Console for a procedure, reference Section 2, “Pre-Operating Instructions”.


Figure 1 – Console Front, Side, and Top Panel


Figure 2 – Console Rear Panel

2. PRE-OPERATING INSTRUCTIONS

Use the steps in this section to prepare the console for a procedure.

2.1. Unpacking and Setup

2.1.1. Console

The Cryoablation Console is supplied in a customized shipping crate. To unpack, the crate must be placed upright on a flat surface.

- For the wooden reusable crate, fold down the crate ramp that is held in place by four (4) external latches. Remove the foam insert at the bottom of the console and unlock the Casters. Carefully roll the console down the ramp.
- For the cardboard disposable crate, cut the two straps using scissors, rotate the four latches on the side of the crate, unwrap the cardboard, remove the foam ramp and secure to base. Remove supporting foams and unlock the Casters. Carefully roll the console down the ramp.

2.1.2. Hose Assemblies

There are two hoses involved in the use of the Console (see Accessories section):

1. The orange 7/8" OD (outer diameter) hose (P/N: 104-1107-001) is a high-pressure hose used to supply N₂ gas to the console and must always be connected during use. It connects to the output valve on top of an available Nitrogen supply tank to the quick disconnect fitting on the back of the Console. This hose has a maximum pressure rating of 3000 PSI.
2. The silver 3" OD hose (P/N: 104-1246-001) is a low-pressure hose used to transfer Liquid Nitrogen refrigerant (LN₂) from an available supply tank into the internal reservoir (Dewar) of the Console via the fill port. This hose has a maximum pressure rating of 175 PSI.

NOTE – All pressure hoses use Diameter Index Safety System (DISS) fittings. The liquid Nitrogen Fill Connector Port is larger than the Pressurized N₂ Gas connection to prevent accidental connection. A country specific connector may be required to adapt hoses to supply tanks outside of the US.

WARNING – PRESSURIZED NITROGEN GAS USED IN CONSOLE – Incorrect Operation may result in damage to the device, ineffective treatment, or harm to the patient and user. Do not attempt to disconnect pressurized Nitrogen Gas line without depressurizing.

2.1.3. Nitrogen Tanks

WARNING – Nitrogen tanks must be properly secured. Gas Nitrogen is under high pressure.

Nitrogen Tanks (Cylinders) are not supplied by Adagio. Always use Ultra High Purity (99.998% pure)/high pressure Nitrogen gas with the Console. The N₂ supply hose has an integrated regulator to reduce the pressure of the console to approximately 800 psi (5.5 N/mm²). No external regulator is necessary on the N₂ tank. Check tank pressures to verify capacity. Install the cylinders and secure.

2.1.4. Handling / Transporting Compressed Gas Cylinder

Gas cylinders must be secured in a vertical position when deployed, stored or transported. Use a chain or strap to fasten cylinders to a solid support, like a wall. Use cylinder trolley or hand truck when moving gas cylinders within the worksite. Always have the cylinder valve cover in place and securely fastened. Close gas cylinder's valve:

- Before moving the cylinder.
- When the work is finished.
- When the cylinder is empty.

Secure in place the valve protection cap until the gas cylinder is ready for use. Valve protection cap should not be used to lift the gas cylinder. The valve protection cap could be damaged and/or become unattached causing the cylinder to drop onto a hard surface, possibly resulting in an explosion.

- Do not use gas cylinders as support or rollers for moving other materials.
- Do not alter numbers and markings stamped on the gas cylinders.
- Mark empty gas cylinders as "Empty".

2.1.5. Storage of Compressed Gas Cylinders

Compressed gas cylinders should be stored in accordance to safety guidelines.

- Compressed gas cylinders are heavy. It is important to handle them with care to prevent dropping, dragging or subjecting them to impact. Before handling or using gas cylinders, understand the properties and hazards of the gas.
- Store cylinders upright and secure the cylinders by the body with straps to a stationary building support or to a cylinder cart to prevent cylinders from tipping or falling over.
- Store cylinders with the valve protection cap fully fastened to the threads.
- Maintain storage areas dry, well ventilated and away from heat sources. A damp storage condition can corrode the gas cylinders.
- Securely connect valve protection cap to gas cylinders when they are not in use.
- Clearly tag/mark empty gas cylinders to indicate their condition.
- Segregate empty gas cylinders from full gas cylinders.
- Store cylinders with their valves in the fully "closed" position.
- Rotate the stock of stored cylinders by using a first-in, first-out (FIFO) inventory system. This is to avoid storing cylinders for an extended period.

2.1.6. Use of Personal Protective Equipment (PPE)

Eyes, hands and feet should be protected.

Wear appropriate PPE when handling compressed gas cylinders. These may include gloves, suitable eye and face protection such as safety glasses/visor and safety shoes.

2.1.7. Cryogen Tank

WARNING – AVOID CONTACT WITH LIQUID NITROGEN CRYO FLUID – Direct contact with Liquid Nitrogen can cause serious burns to exposed skin and eyes. Gloves and goggles should be worn when handling.

Always use Low Pressure (≤ 22 psi / 0.2 N/mm²) Liquid Nitrogen (LN2) to fill the internal reservoir.

2.1.8. Handling Liquid Nitrogen

Contact of liquid nitrogen or any very cold gas with the skin or eyes may cause serious freezing (frostbite) injury.

- Protect hands at all times with safety gloves when working with liquid nitrogen.
- Handle liquid nitrogen carefully and do not dispense from tank rapidly.

The extremely low temperature can freeze human flesh very quickly. When spilled on a surface the liquid tends to cover it completely and intimately, cooling a large area. The gas issuing from the liquid is also extremely cold. Delicate tissue, such as that of the eyes, can be damaged by an exposure to the cold gas which would be too brief to affect the skin of the hands or face.

- Never allow any unprotected part of your body to touch objects cooled by liquid nitrogen.

Such objects may stick fast to the skin and tear the flesh when you attempt to free yourself. Use tongs, preferably with insulated handles, to withdraw objects immersed in the liquid, and handle the object carefully.

- Wear protective clothing.

Protect your eyes with a face shield or safety goggles (safety glasses without side shields do not give adequate protection). Always wear safety gloves when handling anything that is, or may have been, in immediate contact with liquid nitrogen. The gloves should fit loosely, so that they can be thrown off quickly if liquid should splash into them. When handling liquid in open containers, it is advisable to wear high-top shoes. Trousers (which should be cuffless if possible) should be worn outside the shoes.

Any kind of canvas shoes should be avoided because a liquid nitrogen spill can be taken up by the canvas resulting in a far more severe burn, in fact that would occur if the feet were essentially open or bare. It is not recommended to go bare foot or wear canvas when using liquid nitrogen.

- Use only containers designed for low temperature liquids.

Cryogenic containers are specifically designed and made of materials that can withstand the rapid changes and extreme temperature differences encountered in working with liquid nitrogen. Even these special containers should be filled slowly to minimize the internal stresses that occur when any material is cooled. Excessive internal stresses can damage the container.

Do not cover or plug the entrance opening of any liquid nitrogen Dewar. Do not use any stopper or other device that would interfere with venting of gas.

These cryogenic liquid containers are generally designed to operate with little or no internal pressure. Inadequate venting can result in excessive gas pressure which could damage or burst the container. Use only the loose-fitting necktube core supplied or one of the approved accessories for closing the necktube. Check the unit periodically to be sure that venting is not restricted by accumulated ice or frost.

- Use proper transfer equipment.

Use only the supplied Liquid Nitrogen (LN2) hose (P/N: 104-1246-001) when transferring liquid nitrogen into the console Dewar. Be certain both ends of the hose are securely connected before opening valve to initiate transfer.

- Do not overfill containers.

Filling above the bottom of the necktube (or specified maximum level) can result in overflow and spillage of liquid when the necktube core or cover is placed in the opening.

- Liquid nitrogen gas/vapors can cause suffocation without warning. Store and use liquid nitrogen only in a well-ventilated area. As the liquid evaporates, the resulting gas tends to displace the normal air from the area. In closed areas, excessive amounts of nitrogen gas reduce the concentration of oxygen and can result in asphyxiation. Because nitrogen gas is colorless, odorless and tasteless, it cannot be detected by the human senses and will be breathed as if it were air. Breathing an atmosphere that contains less than 19 percent oxygen can cause dizziness and quickly result in unconsciousness and death.

2.1.9. Note

The cloudy vapor that appears when liquid nitrogen is exposed to the air is condensed moisture, not the gas itself. The gas actually causing the condensation and freezing is completely invisible.

Never dispose of liquid nitrogen in confined areas or places where others may enter.

Disposal of liquid nitrogen should be done outdoors in a safe place. Pour the liquid slowly on gravel or bare earth where it can evaporate without causing damage. Do not pour the liquid on the pavement.

- First Aid Notice.

If a person seems to become dizzy or loses consciousness while working with liquid nitrogen, move to a well-ventilated area immediately. If breathing has stopped, apply artificial respiration. If breathing is difficult, give oxygen. Call a physician. Keep warm and at rest.

If exposed to liquid or cold gas, restore tissue to normal body temperature of 98.6°F (37°C) as rapidly as possible, followed by protecting the injured tissue from further damage and infection. Remove or loosen clothing that may constrict blood flow to the frozen area. Call a physician. Rapid warming of the affected part is best achieved by using water at 108°F/42°C. Under no circumstances should the water be over 112°F/44°C, nor should the frozen part be rubbed either before or after rewarming.

Keep the unit clean and dry. Do not store it in wet, dirty areas. Moisture, animal waste, chemicals, strong cleaning agents and other substances which could promote corrosion should be removed promptly. Use water or mild detergent for cleaning and dry the surface thoroughly. Do not use strong alkaline or acid cleaners that could damage the finish and corrode the metal shell.

2.2. Transporting, Positioning and Locking the Console Wheels

Before transporting the Cryoablation Console the following precautions must be observed:

1. Follow the GUI instructions on console shutdown and disconnect the catheter, if connected.
2. Insert the storage plug (104-1247-002) in the catheter connector port.
3. Close the valve on the N2 supply tank.
4. Disconnect the high pressure N2 supply hose & pressure transducer connector from the console.
5. If moving the console only a short distance within the Cath Lab or hospital, it is safe to have LN2 in the console reservoir. Close the fill port valve and use care when moving so as not to slosh the LN2. If the console is to be transported by truck or similar means, the internal Liquid Nitrogen Reservoir must be empty. Allow sufficient time for the Liquid Nitrogen to evaporate (approximately 2-4 days, depending on initial volume).
6. Press and release the On/Off Pushbutton to power down the console.
7. Turn off the Main AC Power Switch on the back panel.
8. Disconnect the AC Power cord from the wall outlet and wrap it around AC power cord bracket on the back panel.
9. Rotate to the rear and fold the touch screen monitor down until it stops.
10. Unlock the locks on the rear casters by moving the caster lock tabs up to the horizontal position.


Figure 3 – Left: Console unlocked; Right: Console locked

NOTE – During transport of the Cryoablation Console, it must never be left unattended on any slope. The casters may not be relied on to hold the equipment in position on a slope.

11. Position the console on a flat horizontal floor in the desired location in the room such that there is direct access to the Main AC Power Switch and power cord.
12. Once the console is in position, push down on the Caster Lock Tab located over the two back wheels to lock the console wheels/casters.
13. Adjust the screen height and angle.

3. OPERATING INSTRUCTIONS - ADAGIO MEDICAL CRYOABLATION CONSOLE (AF)

3.1. Procedure Overview

This overview pertains to ablation procedures using compatible Adagio Medical, Inc. Cryoablation Catheters for treatment of atrial arrhythmias. Detailed console operating instructions begin at Section 3.2 “Pre-cryoablation steps”. Thoroughly understand the Instructions for Use (IFU) for the iCLAS™ catheter for detailed ablation procedure steps.

Pre-cryoablation steps

- Connect power & N2.
- Start-up system.
- Prepare console
- Prepare and insert steerable sheath into the patient.
- Prepare and insert the catheter into the patient.

Cryoablation steps

- Position the catheter.
- Verify catheter positioning.
- Perform cryoablation.

Post-Cryoablation steps

- Remove the catheter from the patient.
- Remove the steerable sheath from the patient.
- Disconnect catheter from console and insert storage plug.
- Shut down console.

3.2. Pre-Cryoablation Steps

1. **Power Connection** – Plug the AC Power Cord into an outlet and move the Main AC Power Switch to the “On” position. This supplies power to the power supplies but does not start the Console.

Note – The AC Power Cord at the rear of the Console is the Emergency Disconnection Device. Position the Console so that there is unobstructed access to it during operation should the need for emergency access arise.


Figure 4 – Console AC Power Connections and Switch

2. **System Start** – The Console is equipped with a Power Button (On/Off) on the rear panel, located on the back panel of the Console. To start the Console, press the Power Button (On/Off) once momentarily and release. It will illuminate within one second.

Note – Pushing it again will turn the Console off.


Figure 5 – On/Off Pushbutton / Power Button

3. Within a few seconds after the On/Off has been pressed, the GUI screen will display the following boot image (*Figure 6*). Windows will prompt you to login. The username is "User". Refer to Adagio Medical Technical Support team for password.

Figure 6 – Boot Screen Image

4. The system will automatically check the communication among its components. While this is in-progress, the following splash screen (Figure 7) will appear. There is no need for user intervention at this time. If the Adagio Medical logo does not appear, click the task view (Figure 8) and select the Adagio Medical GUI application.


Figure 7 – Splash screen image



Figure 8 –How to select Adagio Medical GUI Application image

5. **Program Start** – After the successful communication check which takes a few more seconds, Cryoablation Console Application Main Menu (Figure 9) will appear:

- | | |
|--|--------------------------|
| 1 Add a new patient | 4 Shut down the software |
| 2 Prepare the console for procedure | 5 User settings screen |
| 3 View maintenance screen (password protected) | 6 About Tab |

Figure 9 – Main Menu Screen Image

6. **Main Menu Buttons** – After the Main Menu Screen is displayed, only certain Main Menu Buttons in Figure 9 are active.
7. **Language Setting**
- Language Setting for System User Interface - Press User Settings (gear icon) in the upper right corner of Main Menu (Figure 9), the User Setting screen (Figure 10) will appear. Select the desired language from the dropdown list and press SAVE button to save the selection. Press Back button to return to Main Menu.


Figure 10 – Language Setting for System User Interface Screen Image

- Language Setting for System Keyboard - Once the application launches, slide from the bottom of the screen to bring up task bar.


Figure 11 – Language Setting for System Keyboard Screen Image

3.2.1. “ABOUT TAB”

Open the “ABOUT TAB” to bring up a screen with additional information about software versions of GUI, Main Board, Patient Board, and TCB. NOTE – This information normally is not needed during operation. Click “OK” button to return to the Main Menu screen.

Figure 12 – ABOUT TAB Screen Image

3.3. “PREPARE CONSOLE” FOR A PROCEDURE

Before using the Console, the following steps should be performed in order to load or refill the Liquid N2 Refrigerant into the Console and connect the N2 Gas Supply to provide operating pressure. These steps will be presented in greater details by a series of on-screen animations to demonstrate user actions and ensure safe operations.

- 1. From the Main Menu (Figure 9), select “Prepare Console”. The storage plug must be inserted at this point (Figure 13). Press “CONFIRM” once the storage plug is put in place.


Figure 13 – Insert Storage Plug Screen Image

2. **Connect the N2 Gas Line to the back of the Console.** As shown in *Figure 14* and on-screen animation, connect the pressurized N2 Gas Supply hose to the Console by connecting to the quick-connect port on the back of the Console. Also, connect the Tank Pressure transducer connector into the receptacle (*Figure 14* and *Figure 15*). Select “CONFIRM” upon completion.


Figure 14 – Connect N2 Gas Line To The Console Screen Image

This N2 supply hose must remain connected to the Console to provide operational and functional pressure.

WARNING – PRESSURIZED NITROGEN GAS USED IN CONSOLE – Incorrect Operation may result in damage to the device, ineffective treatment, or harm to the patient and user. Never attempt to disconnect pressurized Nitrogen Gas Line until it is depressurized.

NOTE – If it becomes necessary to disconnect the pressurized Nitrogen Gas Line, it must be first be depressurized. Follow instruction in section 3.6.2 to depressurize the N2 gas line.

3. **Connect the N2 gas tank pressure sensor to the Console** (*Figure 15*) and select “CONFIRM” upon completion. Do not force connect. Be sure pins and holes are aligned before pushing in to connect.


Figure 15 – Connect N2 Gas Tank Pressure Sensor Screen Image

4. The console will automatically run the Prepare Console (Vacuum) test as shown in *Figure 16*.


Figure 16 – Prepare Console Test Screen Image

5. If the console passes the test, *Figure 17* will appear; otherwise, depending on the cause, *Figure 18* or *Figure 19* will be shown. Follow instructions on screen if the Prepare Console Test fails. If the tests pass, press “MAIN MENU” to skip refilling LN2 tank and return to Main Menu.

Figure 17 – Prepare Console Test Passed Screen Image

Figure 18 – Prepare Console Test Failed (Not Enough Vacuum) Screen Image


Figure 19 – Prepare Console Test Failed (Not Enough Pump Speed) Screen Image

6. To refill LN2, from *Figure 17* screen, press “REFILL” to fill liquid nitrogen tank. The console will start purging as needed, to remove any residual moisture from the system (*Figure 20*).


Figure 20 –Console Purge In-progress Screen Image

7. Upon completion of purging, the console LN2 reservoir is ready to be filled. Open the hatch on top of the console, to access the LN2 fill port (*Figure 21*).

Figure 21 – Connect The Liquid N2 Hose To The Console And Fill Its Internal Reservoir Screen Image

WARNING – AVOID CONTACT WITH LIQUID NITROGEN CRYO FLUID – Direct contact with Liquid Nitrogen can cause serious burns to exposed skin and eyes. Safety gloves and goggles should be worn when filling the console LN2 Reservoir.

8. Open the Hatch on top of the Console and open the Valve on the Liquid N2 Fill Port (*Figure 22*).

Fill Valve Closed

Fill Valve Open

Figure 22 – Console Fill Valve

9. Connect the Liquid N2 Hose from the tank by inserting the Fill Tube into the port (*Figure 23*).


Figure 23 – Inserting The Fill Tube Into The Fill Port

10. Pull back on the Fill Tube collar (*Figure 24*) while pushing the Fill Tube into the port. The collar should click when fully seated. After the Fill Tube Collar is locked into the Fill Port, pull up on Fill Tube to check that hose is firmly connected.

Figure 24 – Pulling Back On The Fill Tube Collar To Lock Into The Fill Port

11. In the filling final configuration (*Figure 25*), open the LN2 supply tank valve slowly to initiate the flow of cryogen fluid. Verify no spillage or leaks and gradually open the valve further to continue filling the console reservoir.

Figure 25 – Filling Final Configuration

WARNING – AVOID CONTACT WITH LIQUID NITROGEN CRYO FLUID – Monitor the filling of the Cryogen Fluid and stop filling when it's full. Overfilling the Cryogen Fluid causes it to spill out of the Console onto the floor and potentially contact exposed skin.

12. Once the tank is full, press “CONFIRM” button on *Figure 21* screen. The Refill Complete screen will show (*Figure 26*). Be sure to close the supply tank valve to stop the flow of LN2. Disconnect the fill hose and close the fill port valve. Close the latch cover and then press “MAIN MENU” button on *Figure 26* screen to go back to Main Menu (*Figure 9*).

Figure 26 – Refill Complete screen

3.4. “NEW PATIENT” AND CONNECT CATHETER

The following steps should be performed to connect a Cryoablation Catheter to the Console prior to inserting the Catheter into the patient:

1. From the Main Menu Screen (*Figure 9*), click on the “NEW PATIENT” button to go to Patient Information screen (*Figure 27*).
2. Optionally enter all procedure information into the fields as prompted. Case number is the chronological order of the procedure in a specific date and shall not have any link to patient ID. Upon completion, click the “CONTINUE” button to move on to the next step.

Patient Information

BACK

Case Number *

Physician Name

Center Name

Study Name

Study Date: [] [##] [####]

Comments

* Required field

CONTINUE

Figure 27– Patient Information Screen Image

3. The progress screen indicating vacuum pump status will show next (*Figure 28*).


Figure 28 – Progress Screen Indicating Vacuum Pump Status Screen Image

4. Select the Catheter for use in the procedure as shown in *Figure 29*.


Figure 29 – Catheter Selection Screen Image

5. To confirm your Catheter choice and proceed, click on the “CONFIRM” button.

NOTE – If at any point in the below steps the Catheter fails a test or is not functioning as intended, stop the procedure, disconnect the Catheter, and set the Catheter aside to be reported and returned to Adagio Medical, Inc. Replace defective Catheter with a new one and proceed again starting at main menu.

6. The Console will automatically unlock the Catheter Connection Port.

7. Connect Catheter as displayed on the Screen (Figure 30). If the Catheter fails to connect, click on the “DISCONNECT” button.
- NOTE** – Prior to inserting the Catheter, perform visual inspection per Catheter IFU.


Figure 30 – Catheter Connection Prompt Screen Image

8. After the catheter is successfully connected, the Console will automatically begin to test catheter’s vacuum (Figure 31). If the catheter vacuum test passes, the console will automatically proceed to ask for the catheter to be flushed as in Figure 34 screen.


Figure 31 – Test Catheter Vacuum Screen Image

NOTE – if the catheter vacuum test or connect catheter fails, follow instruction in Figure 32 and Figure 33 to disconnect and replace catheter.


Figure 32 – Catheter Vacuum Test Fail (vacuum and output pressure time out) Screens Image


Figure 33 – Connect Catheter Fail (catheter lock time out) Screens Image

9. To flush the catheter, inject sterile saline as displayed on the Screen (Figure 34). Refer to the Catheter's Instructions for Use for the injection procedure. After the catheter is injected with saline, disconnect the syringe and click on the "PROCEED" button.


Figure 34 – Inject Sterile Saline Screen Image

10. After pressing "PROCEED" in Figure 34, monitor catheter pressure on screen. Follow instruction on animation (Figure 35) to close stopcock and press "PROCEED" to start Catheter self-tests.


Figure 35 – Close Stopcock Screen Image

Note - After saline injection, if the catheter pressure does not meet the pre-set requirement (Figure 36 and Figure 38), the system will ask to open and close the stopcock to adjust catheter pressure (Figure 37) before pressing "PROCEED" to catheter self-tests. There is automated timeout for performing the test and if the test was not performed within 10 minutes, the console will prompt to press disconnect to restart console preparation. Retry saline injection and if it fails, replace the catheter.

Note - This step will allow to relieve residual pressure if the syringe was removed incorrectly before repositioning of the stopcock valve.


Figure 36 – Catheter Pressure Test Failed (Saline Injection Time Out) Screen Image


Figure 37 – Open and Close Stopcock Screen Image

Note - if the Catheter pressure value doesn't meet the acceptable range after the opening and closing stopcock step above, Catheter pressure test failed screen will appear. Press DISCONNECT button and replace the Catheter.


Figure 38 – Catheter Pressure Test Failed (Saline Line Pressure Not Stable) Screen Image

11. The Console will begin to perform a series of Catheter Self-Tests (Figure 39).

NOTE – During these tests the Cryoablation Element (Distal Tip) of the Cryoablation Catheter will briefly become extremely cold (up to -196°C); position it so that it will not cause harm when it does.


Figure 39 – Catheter Self-Test In-progress Screen Image

12. Upon completion of Catheter Self-Tests (Figure 40), user will be prompted with the below message. Click on the “PROCEED” button to enter Full Functional Test mode (Figure 43).


Figure 40 – Catheter Self-Test Complete Screen Image

NOTE – if the catheter self-test fails, follow instruction in Figure 41 and Figure 42 to disconnect and replace catheter.

Figure 41 – Catheter Self-Test Fail (Output Pressure Relief Failed) Screen Image

Figure 42 – Catheter Self-Test Fail (Output Pressure Hold Failed) Screen Image

13. Prior to performing Full Functional Test, immerse the Cryoablation Element (freezing portion) into saline. When ready, in the Full Functional Test screen (Figure 43), click on the “FUNC TEST” button to begin freezing which takes 30 seconds. During the freeze, observe the catheter to ensure that ice forms along the Cryoablation Element and no gas bubbles are released or formed in the ice. In the Full Functional Test In-progress Screen (Figure 44), also observe the Catheter pressure on-screen does not rise above 1 psi. DO NOT USE catheter if bubbles or catheter pressure increase were seen.


Figure 43 – Full Functional Test Screen Image


Figure 44 – Full Functional Test In-progress Screen Image

14. Upon completion of the Full Functional Test, click on the “PROCEED” button to conclude Full Functional Test mode (Figure 45).


Figure 45 – Full Functional Test Complete Screen Image

15. At this time, the Catheter is ready to be inserted into the patient as per the Catheter’s Instructions for Use and the Procedure Ready Screen (Figure 46) will appear.

NOTE – A guide catheter (sheath) must always be used to position the Catheter inside the patient. After the Catheter is positioned, the proximal portion should be clipped to the table or drape to prevent it from accidentally being dislodged.


Figure 46 – Procedure Ready Screen Image

3.5. Cryoablation Steps

3.5.1. Procedure Screen Numerical Values

Terminology	Description
LN ₂ (%)	% of Liquid Nitrogen
N ₂ (%)	% of nitrogen in Dewar
T _D	Distal Temp: Temperature in Catheter Tip (TCD)
T _{IN}	Input Temp: Temperature of NCN flowing into the catheter
T _{OUT}	Output Temp: Temperature of NCN flowing out from the catheter into the console
FLOW	N ₂ gas flow
P _{CAT}	Catheter pressure from a pressure sensor in the catheter handle
V _{AC}	Vacuum in the Catheter

3.5.2. Initiating a Freeze

The Console automatically performs self-checks of critical parameters such as vacuum pressure, gas pressure, refrigerant level, Catheter attachment, etc. before activating the FREEZE button on the Main Menu Screen.

The user should perform the following steps to initiate a Freeze Cycle on the Console:

1. Catheter should be positioned as per Catheter Instructions for Use. After positioning the Catheter, confirm the Procedure Ready screen (Figure 46) is displayed on the monitor.

2. If Cryomapping is desired, press the “CRYOMAPPING” Button on the Procedure Ready screen (Figure 46). The system will go to Cryomapping Screen (Figure 47).


Figure 47 – Cryomapping Screen Image

3. If the user finishes a cryomapping cycle (30 seconds), the system will automatically transition to Freeze screen to continue with a standard freeze cycle as in Figure 46.

The “STOP” Button (Figure 47) allows users to disable the Cryomapping mode prior to timeout (30 seconds). The system will go to the Thaw Screen as in Figure 48.

Once the thaw is complete, the system will go to Thaw Complete screen (Figure 49), press “PROCEED” button to transfer back to Procedure Ready screen (Figure 46) and continue the procedure.


Figure 48 – Catheter Thawing In-progress Screen Image


Figure 49 – Catheter Thawing Complete Screen Image

4. Press the FREEZE Button on the Procedure Ready Screen (Figure 46). At this point the Console will initiate the flow of NCN into the catheter and maintain it for the specified duration. Catheter Freezing Screen (Figure 50) will appear on the monitor. Note that in this screen, the FREEZE button is replaced with a STOP button; if it is pressed again the Freeze Cycle will be terminated prior to the timeout. The Freeze Duration can be also extended by pressing the “ADD CYCLE” button during the Freeze.


Figure 50 – Catheter Freezing Screen Image

5. Carefully observe the Timer and parameters displayed on the Catheter Freezing Screen (Figure 50).
- WARNING – CATHETER ABLATES TISSUE** – The application of cryogenic temperature results in permanent damage to cardiac tissues. Ensure that the cryoablation element is in contact only with the target tissue prior to initiating energy delivery.
- The user can stop the freezing anytime during the cycle by clicking the “STOP” button. The system will transfer to Catheter Thawing screen (Figure 48) and Catheter Thawing Complete screen (Figure 49) once thawing is finished.

6. The catheter pressure (P_{CAT}) will also be monitored by the console software. If P_{CAT} rises to 0.4 psi, the P_{CAT} value will show in orange (Figure 51) at which point the user should pay close attention to P_{CAT} . If P_{CAT} rises to 0.7 psi, the P_{CAT} value will show in red (Figure 52) at which point the user should replace the Catheter. If P_{CAT} rises to 1.0 psi, the Console will eject the Catheter from the Catheter Connection Port. The Storage Plug must be inserted to proceed. Perform visual inspection and Full Functional Test on the Catheter per Catheter IFU.


Figure 51 – P_{CAT} Value In Orange During Freeze Screen Image (Note: P_{CAT} at 0.5 psi is higher than alert limit of 0.4 psi)


Figure 52 – P_{CAT} Value In Red During Freeze Screen Image (Note: P_{CAT} at 0.8 psi is higher than alert limit of 0.7 psi)

7. Upon Freeze completion, the Screen will prompt the user that the catheter is thawing. Once the Catheter Thawing Complete screen (Figure 49) appears, wait until the Thaw is complete (by following the graph on the screen) then you may press “PROCEED” to go back to Procedure Ready Screen (Figure 46); from this screen, you can view procedure summary (Figure 54) and finish the Thaw cycle.

8. The Console can detect when the Catheter has been removed from the patient at any point during the procedure (except during a Freeze cycle). If the Console detects Catheter removal, a dialog box will appear as in (Figure 53) asking the user to select “YES” or “NO”. If “NO” is selected, the dialog box will disappear, and the user can proceed to the next freeze. If “YES” is selected, the Full Functional Test screen (Figure 43) will appear and a Full Functional Test must be performed.


Figure 53 – System Detect Catheter Outside of Patient Screen Image

3.5.3. Terminating Freeze

There are three ways to terminate the Freeze Cycle:

1. The Console will time out the preselected Freeze Cycle Duration and will automatically terminate the Freeze Cycle when that time expires. The Time Indicator Segments will turn blue sequentially to visually indicate the freeze duration; when they all become blue the Freeze will be terminated.
2. After the FREEZE Button has been pressed, it is replaced by the STOP Button. If it's pressed again it will cause the Freeze Cycle to be terminated regardless of the amount of time remaining on the timer.
3. E-Stop Button will automatically terminate Freeze Cycle.

When the Freeze cycle is terminated, continue to carefully observe the Time and parameters displayed on the monitor. Be mindful that if the Catheter Temperature is below 0°C, an ice ball may exist on the Cryoablation Element, thus the Catheter may adhere to Cardiac tissue. Allow sufficient thaw time before removing the Catheter.

WARNING – CATHETER ADHESION – During a Freeze Cycle the Catheter will form an Ice Ball that will adhere to Cardiac Tissue. Do not attempt to manipulate the Catheter during or after cryoablation cycle until sufficient evidence for lack of adhesion to tissue or guide sheath is observed by ensuring the Catheter Distal Temperature (T_D) reaches a minimum of 30°C, as indicated by color change from blue to white, the Thaw cycle as defined by the console software is complete and interference seen on the EGMs during the freeze disappears. Do not rotate, pull, remove, or otherwise manipulate a Catheter that resists movement.

3.6. Post-Cryoablation Steps

3.6.1. Disconnecting a Catheter

- 1. When it becomes necessary to disconnect a catheter, from the Procedure Ready Screen (Figure 46), click “VIEW SUMMARY”. The monitor will display the Procedure Summary Screen (Figure 54).


Figure 54 – Procedure Summary Screen Image

- 2. Press the “END PROCEDURE” button (Figure 54) to disconnect the catheter. The system will confirm if the user would like to disconnect the catheter (Figure 55). Press “DISCONNECT” to continue ending the procedure.


Figure 55 – Confirm Disconnect Catheter Screen Image

3. Console will depressurize the Catheter and automatically unlock the catheter connection (*Figure 56*).
The Catheter can now be removed by gently pulling the Catheter plug out of the connection port. Replace the Catheter with the Storage Plug (if finished with procedures) as shown in *Figure 57*.


Figure 56 – Console Disconnect Catheter Screen Image


Figure 57 – Console Disconnect Catheter – Insert Catheter Plug Screen Image

4. Once the catheter has been removed, replace with the Storage Plug. Navigate to the Main Menu and select “Prepare Console” and repeat section 3.3 steps 1 – 5 and section 3.4.

3.6.2. Console Shutdown

When all procedures are completed, console is to be powered down. Perform the following steps:

1. Confirm to end the console session by pressing “END SESSION” in the below screen *Figure 58*.


Figure 58 – Console Shutdown – Confirmation to End Session Screen Image

2. Shut off the Supply Valve on the External N2 Gas supply tank as prompted in Figure 59 then press “PROCEED”. The console will automatically depressurize the N2 gas line (*Figure 60*).


Figure 59 – Console Shutdown – Turn Off N2 Gas Tank Valve Screen Image

Figure 60 – Console Shutdown – Waiting for Depressurization of N2 Gas Line Screen Image

3. Depressurization of the N2 gas in the system makes noise. When the input pressure reads 0 on the display, it is safe to disconnect the N2 Gas supply hose from the back of console.
4. Disconnect the N2 gas supply hose from the N2 supply tank and store with the console.
5. At this time, the console will return to Main Menu screen (*Figure 9*), press the “SHUT DOWN” button on this Screen.
6. Press and release the On/Off Pushbutton. When the On/Off Pushbutton is darkened, turn off the Main AC Power Switch.

4. OPERATING INSTRUCTIONS - ADAGIO MEDICAL CRYOABLATION CONSOLE (VT)

4.1. Procedure Overview

This overview pertains to ablation procedures using compatible Adagio Medical, Inc. Cryoablation Catheters for treatment of ventricular arrhythmias. Detailed console operating instructions begin at Section 4.2 “Pre-cryoablation steps”. Thoroughly understand the Instructions for Use (IFU) for the vCLAS™ catheter for detailed ablation procedure steps.

Pre-cryoablation steps

- Connect power & N2.
- Start-up system.
- Prepare console.
- Prepare and insert steerable sheath into the patient.
- Prepare and insert the catheter into the patient.

Cryoablation steps

- Position the catheter.
- Verify catheter positioning.
- Perform cryoablation.

Post-Cryoablation steps

- Remove the catheter from the patient.
- Remove the steerable sheath from the patient.
- Disconnect catheter from console and insert storage plug.
- Shut down console.

4.2. Pre-Cryoablation Steps

1. **Power Connection** – Plug the AC Power Cord into an outlet and move the Main AC Power Switch to the “On” position. This supplies power to the power supplies but does not start the Console.

Note – The AC Power Cord at the rear of the Console is the Emergency Disconnection Device. Position the Console so that there is unobstructed access to it during operation should the need for emergency access arise.


Figure 61 – Console AC Power Connections and Switch

2. **System Start** – The Console is equipped with a Power Button (On/Off) on the rear panel, located on the back panel of the Console. To start the Console, press the Power Button (On/Off) once momentarily and release. It will illuminate within one second.

Note – Pushing it again will turn the Console off.


Figure 62 – On/Off Pushbutton / Power Button

3. Within a few seconds after the On/Off has been pressed, the GUI screen will display the following boot image (Figure 63). Windows will prompt you to login. The username is “User”. Refer to Adagio Medical Technical Support team for password.

Figure 63 – Boot Screen Image

4. The system will automatically check the communication among its components. While this is in-progress, the following splash screen (Figure 64) will appear. There is no need for user intervention at this time. If the Adagio Medical logo does not appear, click the task view (Figure 65) and select the Adagio Medical GUI application.


Figure 64 – Splash screen image



Figure 65 –How to select Adagio Medical GUI Application image

5. **Program Start** – After the successful communication check which takes a few more seconds, Cryoablation Console Application Main Menu (Figure 66) will appear:

- | | |
|--|--------------------------|
| 1 Add a new patient | 4 Shut down the software |
| 2 Prepare the console for procedure | 5 User settings screen |
| 3 View maintenance screen (password protected) | 6 About Tab |

Figure 66 – Main Menu Screen Image

6. **Main Menu Buttons** – After the Main Menu Screen is displayed, only certain Main Menu Buttons in Figure 66 are active.

4.2.1. “ABOUT TAB”

Open the “ABOUT TAB” to bring up a screen with additional information about software versions of GUI, Main Board, Patient Board, and TCB. NOTE – This information normally is not needed during operation. Click “OK” button to return to the Main Menu screen.

Figure 67 – ABOUT TAB Screen Image

4.3. “PREPARE CONSOLE” for a Procedure

Before using the Console, the following steps should be performed in order to load or refill the Liquid N2 Refrigerant into the Console and connect the N2 Gas Supply to provide operating pressure. These steps will be presented in greater details by a series of on-screen animations to demonstrate user actions and ensure safe operations.

1. From the Main Menu (Figure 66), select “Prepare Console”. The storage plug must be inserted at this point (Figure 68). Press “CONFIRM” once the storage plug is put in place.


Figure 68 – Insert Storage Plug Screen Image

2. **Connect the N2 Gas Line to the back of the Console.** As shown in Figure 69 and on-screen animation, connect the pressurized N2 Gas Supply hose to the Console by connecting to the quick-connect port on the back of the Console. Also, connect the Tank Pressure transducer connector into the receptacle (Figure 69 and Figure 70). Select “CONFIRM” upon completion.


Figure 69 – Connect N2 Gas Line To The Console Screen Image

This N2 supply hose must remain connected to the Console to provide operational and functional pressure.

WARNING – PRESSURIZED NITROGEN GAS USED IN CONSOLE – Incorrect Operation may result in damage to the device, ineffective treatment, or harm to the patient and user. Never attempt to disconnect pressurized Nitrogen Gas Line until it is depressurized.

NOTE – If it becomes necessary to disconnect the pressurized Nitrogen Gas Line, it must be first be depressurized. Follow instruction in section 3.6.2 to depressurize the N2 gas line.

3. **Connect the N2 gas tank pressure sensor to the Console** (Figure 70) and select “CONFIRM” upon completion. Do not force connect. Be sure pins and holes are aligned before pushing in to connect.


Figure 70 – Connect N2 Gas Tank Pressure Sensor Screen Image

4. The console will automatically run the Prepare Console (Vacuum) test as shown in Figure 71.


Figure 71 – Prepare Console Test Screen Image

5. If the console passes the test, *Figure 72* will appear; otherwise, depending on the cause, *Figure 73* or *Figure 74* will be shown. Follow instruction on screen if the Prepare Console Test fails. If the tests pass, press “MAIN MENU” to skip refilling LN2 tank and return to Main Menu.

Figure 72 – Prepare Console Test Passed Screen Image


Figure 73 – Prepare Console Test Failed (Not Enough Vacuum) Screen Image

Figure 74 – Prepare Console Test Failed (Not Enough Pump Speed) Screen Image

6. To refill LN2, from Figure 72 screen, press “REFILL” to fill liquid nitrogen tank. The console will start purging as needed, to remove any warm LN2 from the reservoir (Figure 75).


Figure 75 –Console Purge In-progress Screen Image

7. Upon completion of purging, the console LN2 reservoir is ready to be filled. Open the hatch on top of the console, to access the LN2 fill port (Figure 76).

Figure 76 – Connect The Liquid N2 Hose To The Console And Fill Its Internal Reservoir Screen Image

WARNING – AVOID CONTACT WITH LIQUID NITROGEN CRYO FLUID – Direct contact with Liquid Nitrogen can cause serious burns to exposed skin and eyes. Safety gloves and goggles should be worn when filling the console LN2 Reservoir.

8. Open the Hatch on top of the Console and open the Valve on the Liquid N2 Fill Port (Figure 77).

Figure 77 – Console Fill Valve

9. Connect the Liquid N2 Hose from the tank by inserting the Fill Tube into the port (*Figure 78*).


Figure 78 – Inserting The Fill Tube Into The Fill Port

10. Pull back on the Fill Tube collar (*Figure 79*) while pushing the Fill Tube into the port. The collar should click when fully seated. After the Fill Tube Collar is locked into the Fill Port, pull up on Fill Tube to check that hose is firmly connected.

Figure 79 – Pulling Back On The Fill Tube Collar To Lock Into The Fill Port

11. In the filling final configuration (*Figure 80*), open the LN2 supply tank valve slowly to initiate the flow of cryogen fluid. Verify no spillage or leaks and gradually open the valve further to continue filling the console reservoir.

Figure 80 – Filling Final Configuration

WARNING – AVOID CONTACT WITH LIQUID NITROGEN CRYO FLUID – Monitor the filling of the Cryogen Fluid and stop filling when it's full. Overfilling the Cryogen Fluid causes it to spill out of the Console onto the floor and potentially contact exposed skin.

12. Once the tank is full, press “CONFIRM” button on *Figure 76* screen. The Refill Complete screen will show (*Figure 81*). Be sure to close the supply tank valve to stop the flow of LN2. Disconnect the fill hose and close the fill port valve. Close the latch cover and then press “MAIN MENU” button on *Figure 81* screen to go back to Main Menu (*Figure 66*).

Figure 81 – Refill Complete screen

4.4. “NEW PATIENT” and Connect Catheter

The following steps should be performed to connect a Cryoablation Catheter to the Console prior to inserting the Catheter into the patient:

- 1. From the Main Menu Screen (Figure 66), click on the “NEW PATIENT” button to go to Patient Information screen (Figure 82).
- 2. Optionally enter all procedure information into the fields as prompted. Case number is the chronological order of the procedure in a specific date and shall not have any link to patient ID. Upon completion, click the “CONTINUE” button to move on to the next step.


Figure 82 – Patient Information Screen Image

- 3. Select the Catheter for use in the procedure as shown in Figure 83.


Figure 83 – Catheter Selection Screen Image

- 4. To confirm your Catheter choice and proceed, click on the “CONFIRM” button.
NOTE – If at any point in the below steps the Catheter fails a test or is not functioning as intended, stop the procedure, disconnect the Catheter, and set the Catheter aside to be reported and returned to Adagio Medical, Inc. Replace defective Catheter with a new one and proceed again starting at main menu.
- 5. The Console will automatically unlock the Catheter Connection Port.

6. Connect Catheter as displayed on the Screen (Figure 84). If the Catheter fails to connect, click on the “DISCONNECT” button.


Figure 84 – Catheter Connection Prompt Screen Image

7. After the catheter is successfully connected, the Console will automatically begin to test catheter’s vacuum (Figure 85). If the catheter vacuum test passes, the console will automatically proceed to ask for the catheter self-test as in Figure 88 screen.


Figure 85 – Test Catheter Vacuum Screen Image

NOTE – if the catheter vacuum test or connect catheter fails, follow instruction in Figure 86 and Figure 87 to disconnect and replace catheter.

Figure 86 – Catheter Vacuum Test Fail (vacuum and output pressure time out) Screen Image


Figure 87 – Connect Catheter Fail (catheter lock time out) Screen Image

8. After the catheter connection test is successful, the Console will begin to perform a series of Catheter Self-Tests (*Figure 88*).
NOTE – During these tests the Cryoablation Element (Distal Tip) of the Cryoablation Catheter will briefly become extremely cold (up to -196°C); position it so that it will not cause harm when it does.


Figure 88 – Catheter Self-Test In-progress Screen Image

9. Upon completion of Catheter Self-Tests (*Figure 89*), user will be prompted with the below message. Click on the “PROCEED” button to enter Full Functional Test mode (*Figure 92*).

Figure 89 – Catheter Self-Test Complete Screen Image

NOTE – if the catheter self-test fails, follow instruction in Figure 90 and Figure 91 to disconnect and replace catheter.


Figure 90 – Catheter Self-Test Fail (Output Pressure Hold Failed) Screen Image

Figure 91 – Catheter Self-Test Fail (Output Pressure Relief Failed) Screen Image

10. Prior to performing Full Functional Test, immerse the Cryoablation Element (freezing portion) into sterile saline. When ready, in the Full Functional Test screen (Figure 92), click on the “FUNC TEST” button to begin freezing which takes 30 seconds. During the freeze, observe the catheter to ensure that ice forms along the Cryoablation Element and no gas bubbles are released or formed in the ice. Also observe that the Catheter proximal temperature on-screen remains above -10°C. The lack of ice formation, observed bubbles, or Catheter Proximal Temperature (T_p) below -10°C represent the failure of Full Functional Test. The catheter MUST be replaced.

NOTE – During these tests the Cryoablation Element (Distal Tip) of the Cryoablation Catheter will briefly become extremely cold (up to -190°C); keep the distal tip submerged in sterile saline at all times during this test.


Figure 92 – Full Functional Test Screen Image


Figure 93 – Full Functional Test In-progress Screen Image

11. Upon completion of the Full Functional Test (Figure 94), click on the “PROCEED” button to conclude Full Functional Test mode.


Figure 94 – Full Functional Test Complete Screen Image

12. At this time, the Catheter is ready to be inserted into the patient as per the Catheter’s Instructions for Use and the Procedure Ready Screen (Figure 95) will appear.

NOTE – A guide catheter (sheath) must always be used to position the Catheter inside the patient. After the Catheter is positioned, the proximal portion should be clipped to the table or drape to prevent it from accidentally being dislodged.


Figure 95 – Procedure Ready Screen Image

4.5. Cryoablation Steps

4.5.1. Procedure Screen Numerical Values

Terminology	Description
LN ₂ (%)	% of Liquid Nitrogen
N ₂ (%)	% of nitrogen in Dewar
T _D	Distal Temp: Temperature in Catheter Tip (TCD)
T _P	Proximal Temp: Temperature in Catheter Proximal (TCP)
T _{IN}	Input Temp: Temperature of NCN flowing into the catheter
T _{OUT}	Output Temp: Temperature of NCN flowing out from the catheter into the console
FLOW	N2 gas flow
V _{AC}	Vacuum in the Catheter

4.5.2. Initiating a Freeze

The Console automatically performs self-checks of critical parameters such as vacuum pressure, gas pressure, refrigerant level, Catheter attachment, etc. before activating the FREEZE button on the Main Menu Screen.

The user should perform the following steps to initiate a Freeze Cycle on the Console:

1. Catheter should be positioned as per Catheter Instructions for Use. After positioning the Catheter, confirm the Procedure Ready screen (Figure 95) is displayed on the monitor.
2. If Cryomapping is desired, press the "CRYOMAPPING" Button on the Procedure Ready screen (Figure 46). The system will go to Cryomapping Screen (Figure 96).


Figure 96 – Cryomapping Screen Image

3. If the user finishes a cryomapping cycle (30 seconds), the system will automatically transition to Freeze screen to continue with a standard freeze cycle as in Figure 95.

The "STOP" Button (Figure 96) allows users to disable the Cryomapping mode prior to timeout (30 seconds). The system will go to the Thaw Screen as in Figure 97.

Once the thaw is complete, the system will go to Thaw Complete screen (Figure 98), press "PROCEED" button to transfer back to Procedure Ready screen (Figure 95).


Figure 97 – Thaw Screen In-progress Image


Figure 98 – Thaw Complete Screen Image

4. Press the FREEZE Button on the Procedure Ready Screen (Figure 95). At this point the Console will initiate the flow of NCN into the catheter and maintain it for the specified duration. Catheter Freezing Screen (Figure 99) will appear on the monitor. Note that in this screen, the FREEZE button is replaced with a STOP button; if it is pressed again the Freeze Cycle will be terminated prior to the timeout. The Freeze Duration can also be extended by pressing the “ADD CYCLE” Button during the Freeze.


Figure 99 – Catheter Freezing Screen Image

5. Carefully observe the Timer and parameters displayed on the Catheter Freezing Screen (Figure 99).
- WARNING – CATHETER ABLATES TISSUE** – When a Freeze Cycle is in progress, the Cryoablation Element (freezing zone) may reach NCN temperature (-196°C) and will ablate tissue. Be certain the Catheter is in the correct position in the heart before initiating the Freeze Cycle.
- The default duration of the freeze cycle is 30 seconds. Press “ADD CYCLE” to increase freeze duration.
 - CONTINUOUSLY monitor the catheter’s proximal temperature on the Console display during cryogen delivery. A decrease of the Catheter Proximal Temperature (less than -10°C) indicates a potential internal leak. IMMEDIATELY DISCONNECT and REMOVE the Catheter from patient.
 - The user can stop freezing anytime during the cycle by clicking the “STOP” button. The system will transfer to Catheter Thawing screen (Figure 97) and Catheter Thawing Complete screen (Figure 98) once thawing is finished.
6. Upon Freeze completion, the Screen will prompt the user that the catheter is thawing (Figure 100). Once the Catheter Thawing Complete screen (Figure 101) appears, wait until the Thaw is complete (by following the graph on the screen) then you may press “PROCEED” to go back to Procedure Ready Screen (Figure 95). From this screen, you can view procedure summary (Figure 103).


Figure 100 – Catheter Thawing In-progress Screen Image


Figure 101 – Catheter Thawing Complete Screen Image

7. The Console can detect when the Catheter has been removed from the patient at any point during the procedure (except during a Freeze cycle). If the Console detects Catheter removal, a dialog box will appear as in (Figure 103) asking the user to select “YES” or “NO”. If “NO” is selected, the dialog box will disappear, and the user can proceed to the next freeze. If “YES” is selected, the Full Functional Test screen (Figure 92) will appear and a Full Functional Test must be performed.


Figure 102 – System Detect Catheter Outside of Patient Screen Image

4.5.3. Terminating Freeze

There are three ways to terminate the Freeze Cycle:

1. The Console will time out the preselected Freeze Cycle Duration and will automatically terminate the Freeze Cycle when that time expires. The Time Indicator Segments will turn blue sequentially to visually indicate the freeze duration; when they all become blue the Freeze will be terminated.
2. After the FREEZE Button has been pressed, it is replaced by the STOP Button. If it's pressed again it will cause the Freeze Cycle to be terminated regardless of the amount of time remaining on the timer.
3. Quick Stop Button will automatically terminate Freeze Cycle.

When the Freeze cycle is terminated, continue to carefully observe the Time and parameters displayed on the monitor. Be mindful that if the Catheter Temperature is below 0°C, an ice ball may exist on the Cryoablation Element, thus the Catheter may adhere to cardiac tissue. Allow sufficient thaw time before removing the Catheter.

WARNING – CATHETER ADHESION – During a Freeze Cycle the Catheter will form an Ice Ball that will adhere to Cardiac Tissue and may adhere to the guide sheath. Do not attempt to manipulate the Catheter during or after cryoablation cycle until sufficient evidence for lack of adhesion to tissue or guide sheath is observed by ensuring the Catheter Distal Temperature (T_D) reaches a minimum of 30°C, as indicated by color change from blue to white, the Thaw cycle as defined by the console software is complete, and interference seen on the EGMs during the freeze disappears. Do not rotate, pull, deflect, or otherwise manipulate a Catheter that resists movement.

4.6. Post-Cryoablation Steps

4.6.1. Disconnecting a Catheter

- 1. When it becomes necessary to disconnect a catheter, from the Procedure Ready Screen (Figure 95), click “VIEW SUMMARY”. The monitor will display the Procedure Summary Screen (Figure 103).


Figure 103 – Procedure Summary Screen Image

- 2. Press the “END PROCEDURE” button (Figure 103) to disconnect the catheter. The system will confirm if the user would like to disconnect the catheter (Figure 104). Press “DISCONNECT” to continue ending the procedure.


Figure 104 – Confirm Disconnect Catheter Screen Image

3. Console will depressurize the Catheter and automatically unlock the catheter connection.

The Catheter can now be removed by gently pulling the Catheter plug out of the connection port (*Figure 105*).

The Catheter can be replaced with another catheter or with the Storage Plug (if finished with procedures) as shown in *Figure 106*.

Figure 105 – Console Disconnect Catheter Screen Image


Figure 106 – Console Disconnect Catheter – Insert Catheter Plug Screen Image

4.6.2. Console Shutdown

When all procedures are completed, console is to be powered down. Perform the following steps:

1. Confirm to end the console session by pressing “END SESSION” in the below screen (*Figure 107*).


Figure 107 – Console Shutdown – Confirmation to End Session Screen Image

2. Shut off the Supply Valve on the External N2 Gas supply tank as prompted in *Figure 108* then press “PROCEED”. The console will automatically depressurize the N2 gas line (*Figure 109*).


Figure 108 – Console Shutdown – Turn Off N2 Gas Tank Valve Screen Image

Figure 109 – Console Shutdown – Waiting for Depressurization of N2 Gas Line Screen Image

3. Depressurization of the N2 gas in the system makes noise. When the input pressure reads 0 on the display, it is safe to disconnect the N2 Gas supply hose from the back of console.
4. Disconnect the N2 gas supply hose from the N2 supply tank and store with the console.
5. At this time, the console will return to Main Menu screen (*Figure 66*), press the “SHUT DOWN” button on this Screen.
6. Press and release the On/Off Pushbutton. When the On/Off Pushbutton is darkened, turn off the Main AC Power Switch.

5. SPECIFICATIONS

5.1. General Mechanical Specifications¹

	Height inches (cm)	Width inches (cm)	Depth inches (cm)	Weight lbs. (kg)
Console Including Monitor	56 (142.2)	24 (61.0)	25 (63.5)	230 (105)

Table 2 – Console Dimensions

5.2. Electrical Specifications

AC Power-----120 VAC 50-60 Hz 3.6 Amps
 220-240 VAC 50-60 Hz 2.0 – 1.8 Amps

5.3. Operational Specifications and Environment

Temperature-----10°C to 30°C (50°F to 86°F)
 Humidity-----10% to 75% non-condensing
 Pressure-----70 kPa to 106 kPa
 Altitude-----400 to 3000 m (-1312' to 9842')

5.4. Transport Environment

Temperature-----18°C to 38°C
 Humidity-----max 85%
 Pressure-----70 kPa to 106 kPa

5.5. Technical Data

Manufacturer	Adagio Medical Inc. Laguna Hills, CA 92653
Model	Cryoablation System
Classification IEC 60601-1	-Class I -Applied part, Type CF (defibrillation proof)
Mode	Continuous Operation
Supply Voltage and Power Consumption	120V~50-60Hz 3.6A 220-240V~50-60Hz 2.0-1.8A
Fuse / Rating	Two (one per AC leg) Slow-Blow, miniature 5×20 mm, 3.15A 250V (EU), 6.3A 125V (US and Canada)
Console	Equipment not suitable for use in presence of a FLAMABLE ANESTHETIC MIXTURE WITH AIR or NITROUS OXIDE
Pressure	N2 Operating Gas: 3000 PSI Tank 750 PSI First (Hose) Regulator +/- 10%
Flow	Nominal: 60 Standard Liters per Minute. Flow is automatically controlled by the console and cannot be changed by the user.
Weight	230 lbs.
Console Dimensions (H X W X D)	56" × 24" × 25"

Table 3 – Technical Data

¹ Dimensions and weights are approximate and do not include any accessory equipment. Height dimensions are given such that when added the approximate total height is given.

5.6. Console Configuration Kit

Item#	Part Name	Quantity	Part Number
1	Operator's Manual	1	108-0063-001 (EU) 108-0063-002 (US) 108-0063-003 (Canada)
2	Power Cord	1	102-0236-001 (EU) 102-0236-002 (US and Canada)
3	Fuse	1	102-0060-001 (US and Canada) 102-0060-002 (EU)
4	High Pressure Nitrogen Gas (N2) Hose	1	104-1107-001
5	Liquid Nitrogen (LN2) Fill Hose	1	104-1246-001
6	Console Storage Plug	1	104-1247-002
7	Accessories Case	1	101-0762-001
8	2 Adjustable Wrenches (LN2 Hose to Tank)	1	105-0347-001
9	Safety Gloves for LN2 Fill	1	114-0037-001
10	Hex Wrench Set	1	105-0348-001
11	Phillips Screwdriver Set	1	105-0349-001
12	Console Cover	1	107-0064-001
13	Nitrogen Bottle Connection	1	Installed (US and Canada) 101-0771-001 (Netherlands) 101-0772-001 (Germany) 101-0773-001 (France) 101-0815-001 (Belgium)

Table 4 – Accessories

6. MANUFACTURER'S DECLARATION

6.1. Manufacturer's declaration – Electromagnetic Emission (Table 2, IEC 60601-1-2:2014, 4.1 Edition)

The product is suitable for use in a specific electromagnetic environment. The customer and/or the user of the product should assure that it is used in an electromagnetic environment as described below.

Emissions Test	Compliance	Electromagnetic Environment
Radiated RF emission CISPR 11	Complies	The product is limited to use in Professional Healthcare Facility environment.
Conducted RF emission CISPR 11	Complies	
Harmonic emissions IEC 61000-3-2	Complies	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Complies	

Table 5 – Electromagnetic Emission

6.2. Manufacturer's declaration – Electromagnetic Immunity I (Table 4, 5, and 7 IEC 60601-1-2:2014, 4.1 Edition)

The product is suitable for use in a specific electromagnetic environment. The user of the product should assure it is used in an electromagnetic environment described below.

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment
Electrostatic discharge (ESD) IEC 61000-4-2	± 8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air	Complies	The product is limited to use in Professional Healthcare Facility hospital environment.
Electrical fast transient/bursts IEC 61000-4-4	± 2 kV and 100 kHz repetition frequency for power supply lines ± 1 kV and 100 kHz repetition frequency for input/output lines	Complies	Mains power quality should be that of a typical hospital environment
Surge IEC 61000-4-5	± 0.5 kV, ± 1 kV, ± 2 kV (Line-to-ground) ± 0.5 kV, ± 2 kV (Line-to-line)	Complies	Mains power quality should be that of a typical hospital environment
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	Voltage Dips: < 0 % U_T for 0.5 cycle at 0°, 45°, 90°, 135°, 180°, 225°, 210° and 315°. < 0 % U_T for 1 cycle and 70 % U_T for 25/30 cycles single phase at 0° Voltage Interruptions: 0 % U_T ; 250/300 cycle	Complies	Mains power quality should be that of a typical hospital environment. If the user of the product requires continued operation during power mains interruptions, it is recommended that the product be powered from an uninterruptible power supply.
RATED Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m, 50/60 Hz	Complies	Power frequency magnetic fields should be at levels characteristic of a typical location in a hospital environment.

Note: U_T is the mains (AC) voltage before apply test levels

Table 6 – Electromagnetic Immunity

6.3. Manufacturer's declaration – Electromagnetic Immunity II (Table 5, 7, and 8, IEC 60601-1-2:2014, 4.1 Edition)

The product is suitable for use in a specific electromagnetic environment. The user of the product should assure that it is used in an electromagnetic environment as described below.

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment
Conducted Disturbances Induced by RF Fields IEC 61000-4-6	3 V 150 kHz to 80 MHz 6 V in ISM bands between 150 kHz to 80 MHz 80 % AM at 1 kHz	Complies	Portable and mobile RF communications equipment should be used no closer than 30 cm (12 in) to any part of the product, including cables, and the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.
Radiated RF EM IEC 61000-4-3	3 V/m 80 MHz to 2.7 GHz 80 % AM at 1 kHz	Complies	Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey ^a , should be less than the compliance level in each frequency range ^b . Interference may occur in the vicinity of equipment marked with the symbol. 

Note 1: At 80 MHz and 800 MHz, the higher frequency range applies.

Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, people and animals.

^a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered, if the measured field strength in the location in which the product is used exceeds the applicable RF compliance level above, the product should be observed, additional measures may be necessary, such as reorienting or relocating the product.

^b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

Table 7 – Electromagnetic Immunity II

6.4. Manufacturer's declaration – Recommended Separation Distances Between Portable and Mobile RF-Communications Equipment and the Product (Table 9, IEC 60601-1-2:2014, 4.1 Edition)

The product is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the product can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the product – according to output power and frequency of the communications equipment – as recommended in the following table.

Immunity Test	IEC 60601 Test Level	Compliance Level	Separation Distance in Meters
Proximity Fields from RF Wireless Communication Equipment IEC 61000-4-3	See Subclause 8.10 in 60601-1-2:2014, Table 9	Complies	See Subclause 8.10 in 60601-1-2:2014, Table 9

Table 8 – Separation Distance

7. MAINTENANCE AND SERVICING

WARNING – HIGH VOLTAGES – Dangerous voltages may be contained on circuitry internal to the unit. Maintenance should be performed by qualified personnel only.

CAUTION – The Console electronics contain static sensitive components that can be damaged by improper handling. Use approved grounding techniques for work areas and service personnel.

7.1. Preventive Maintenance

Adagio Medical, Inc. recommends that the Console and its Accessories undergo preventative maintenance as indicated below. To schedule this, contact your Adagio Medical, Inc. service representative.

7.1.1. Annual Maintenance

- Check the refrigerant tank gauge reading annually for proper function.
- Check inline gas purifier filter.
- Clean screens on all vents/fans.
- Preventative Maintenance Checks.
- Check Lithium Ion Batteries (P/N: 102-0206-001):

The Console utilizes four (4) rechargeable lithium ion batteries. The lithium ion batteries power the computer and monitor in the event of a power outage. When the Cryoablation Console is not plugged into a wall socket the battery has an estimated life of three years.

Batteries are only manufacturer serviceable and will be disposed of per local government regulations.

NOTE – Service work should be carried out by Adagio Medical, Inc. service personnel only!

7.1.2. 5-year Service Maintenance

7.1.2.1 Diaphragm Pump (MVP-010-3 DC) (P/N: 102-0231-001):

The valves and diaphragm are wear parts. If the rated ultimate vacuum is no longer achieved, return the Console to Adagio Medical, Inc. for valve cleaning and inspection.

In case of normal wear, the lifetime of the diaphragms and valves is 17500 operating hours at normal rotational speed.

Frequency	Servicing
Every 5 years (or 17500 hours) whichever comes first.	Maintenance of Diaphragm Pump

Table 9 – Diaphragm Servicing Frequency

NOTE – Service work should be carried out by Adagio Medical, Inc. service personnel only!

7.1.2.2 Turbopump (HIPACE 10) (P/N: 102-0134-001):

Perform preventative maintenance :

Frequency	Servicing
Every 5 years	Maintenance of Turbo Pump

Table 10 – Turbopump Servicing Frequency

NOTE – Service work should be carried out by Adagio Medical, Inc. service personnel only!

7.2. Software Updates

The system has been designed to not require any updates to the off-the-shelf software components. Any updates to proprietary software will be performed by Adagio trained personnel only.

7.3. Refrigerant Shipping

To ship the refrigerant to other sites, you must adhere to all local and hospital transportation laws for dangerous materials.

7.4. Safe Disposal

The Console including batteries shall be returned to Adagio Medical, Inc. for safe disposal as per local government regulations.

8. INFORMATIONAL MESSAGES, ERRORS, AND TROUBLESHOOTING

8.1. Informational Messages

Informational messages will appear in blinking yellow text at the bottom of the current screen. Up to six (6) messages can be displayed at the same time. The below table presents a list of messages and troubleshooting.

Message Text	Troubleshooting
CATHETER VACUUM ABOVE THRESHOLD	Replace catheter. Note: if this happens during a freeze, the system will display Safe Abort Pop-up as in <i>Table 13</i> and go to safe abort state.
INPUT PRESSURE NOT IN RANGE	Check console connection and gas tank pressure. Replace N2 gas tank if it is empty. Reboot Console. If the problem persists, contact technical support.
DISTAL TEMPERATURE NOT IN RANGE	Replace or reconnect catheter.
PROXIMAL TEMPERATURE NOT IN RANGE	Replace or reconnect catheter. Proceed with caution if reconnecting catheter does not fix the issue.
CATHETER PRESSURE OUT OF RANGE	Replace catheter (iCLAS and/or iCLAS CoreTx Catheter) Note: if this happens during a freeze, the system will display Safe Abort Pop-up as in <i>Table 13</i> and go to safe abort state. The Console will immediately eject the Catheter. Replace it with the Storage plug.
CATHETER PROXIMAL TOO COLD	Replace catheter (vCLAS Catheter) Note: if this happens during a freeze, the system will display Safe Abort Pop-up as in <i>Table 13</i> and go to safe abort state. The Console will immediately eject the Catheter. Replace it with the Storage plug.
OUTPUT PRESSURE NOT IN RANGE	Check console connection and tank pressure. Reboot Console. If the problem persists, contact technical support.
VALVE PRESSURE NOT IN RANGE	Check console connection and tank pressure. Reboot Console. If the problem persists, contact technical support.
TANK PRESSURE NOT IN RANGE	Ensure Tank Pressure sensor is connected, tank is full and valve on tank is open.
LN2 TANK LOW	Refill tank. Follow instruction in section 3.3 to fill LN2 tank (Dewar).
CONSOLE VACUUM NOT IN RANGE	Reboot Console. If the problem persists, contact technical support.
INPUT TEMPERATURE NOT IN RANGE	Reboot Console. If the problem persists, contact technical support.
OUTPUT TEMPERATURE NOT IN RANGE	Reboot Console. If the problem persists, contact technical support.
EXHAUST TEMPERATURE NOT VALID	Proceed with caution. Contact technical support after the procedure.
EXHAUST TEMPERATURE LOW	Exhaust fan automatically turns on to warm the system. Proceed with caution. Contact technical support after the procedure for assistance if the fan does not turn on.
DEWAR SENSOR DISCONNECTED	Proceed with caution.
RETURN TEMPERATURE WARM	Check the Dewar level on the monitor to make sure the tank is not empty. If the tank is not empty, contact technical support. If during the procedure, proceed with caution. Contact technical support after the procedure for assistance.
FLOW OUT OF RANGE	If during the procedure, proceed with caution. Contact technical support after the procedure for assistance.
CATHETER LOCK DISENGAGED	Disconnect catheter and reinsert fully.
OUTPUT PRESSURE TOO HIGH TO UNLOCK	Reboot Console. If the problem persists, contact technical support.
LN2 TANK WARM	Refill Liquid N2. Follow instruction in Section 3.3 to fill liquid nitrogen tank.
LIQUID N2 LATCH OPEN	Close Liquid Nitrogen Latch. Note: if this happens during a freeze, the system will display informational message. If the latch is still not closed, the console will not proceed.
LN2 TANK FULL	During console setup this message appears when the LN2 tank (Dewar) is filled to 100%, stop filling the Dewar.

Message Text	Troubleshooting
EXHAUST TEMPERATURE COLD	Exhaust fan automatically turns on to warm the system. Proceed with caution. Contact technical support after the procedure for assistance if the fan does not turn on.
CATHETER DISCONNECTED	Catheter is disconnected. If this happens during catheter tests or procedure, the console will not proceed.
PB SUPPLIES OUT OF RANGE	Reboot Console. If the problem persists, contact technical support.
CATHETER VACUUM GAUGE NOT CONFIGURED	Reboot Console. If the problem persists, contact technical support.
CONSOLE VACUUM GAUGE NOT CONFIGURED	Reboot Console. If the problem persists, contact technical support.

Table 11 – List of Informational Messages

8.2. Error Pop-ups

When errors happen, an error pop-up will show with only “Shutdown” as the option. Error pop-up means the device is not safe to operate and needs to be shutdown.

An example error pop-up with a “Power On Self Test Failure” is displayed below:


Figure 110 – An Example of Error Pop-up

Other error conditions are listed in the following table:

Error Title	Error Code	Description
Console Error Detected	E-000	Console generic failure.
Power on Self-Test Failure	E-001	Power on Self-Test Failure. Error code will be different depending on the hardware readings and reason for the failure.
Console Error Detected	E-002	Port connections failed: Main and/or Pump
E-Stop button pressed	E-003	E-Stop button is pressed.
Console Error Detected	E-004	Supply 24V is out of range.
Console Error Detected	E-005	Supply 5V is out of range.
Console Error Detected	E-006	Supply 3.3V is out of range.
Console Error Detected	E-007	Supply -5V is out of range.
Console Error Detected	E-008	TC voltage reference is out of range.
Console Error Detected	E-009	TC level board is not responding.
Console Error Detected	E-010	Catheter vacuum gauge is not responding.
Console Error Detected	E-011	Console vacuum gauge is not responding.
Console Error Detected	E-012	Patient Board is not responding.
Console Error Detected	E-013	Main board is not responding.
Console Error Detected	E-014	ADC timed out.
Console Error Detected	E-015	ADC channel is out of range.
Console Error Detected	E-016	Timer config failed.
Console Error Detected	E-017	Input pressure relief timed out.
Console Error Detected	E-018	Console is cold, vacuum not low.

Error Title	Error Code	Description
Console Error Detected	E-019	Main board RTD chip failed.
Console Error Detected	E-020	Main board HAL initialization error.
Console Error Detected	E-021	TC level board initialization failed.
Console Error Detected	E-022	Patient board initialization failed.
Console Error Detected	E-023	Main board watchdog error.
Console Error Detected	E-024	Output pressure is too high to unlock.


Table 12 – List of Errors Pop-ups

8.3. Safe Abort Screens


If any problem is detected with the catheter, safe abort dialog pop-ups indicating there's a problem with the catheter and need to be disconnected. Pressing 'Disconnect' will transition to Disconnect screen.

An example safe abort pop-up with a "Catheter pressure out of range." reason is displayed as below:

For iCLAS or iCLAS CoreTx Catheter

For vCLAS Catheter

*Figure 111– Examples of Error Pop-up*

Possible reasons for safe abort pop-ups are listed in the following table:

Title	Action	Reason
Problem detected with the catheter	Replace Catheter	Input pressure not in range.
Problem detected with the catheter	Replace Catheter	Catheter lock disengaged.
Problem detected with the catheter	Replace Catheter	Catheter disconnected.
Problem detected with the catheter	Replace Catheter	Catheter vacuum above threshold.
Problem detected with the catheter	Replace Catheter	Catheter pressure out of range (iCLAS and iCLAS CoreTx Catheters only).
Problem detected with the catheter	Replace Catheter	Catheter proximal too cold (vCLAS Catheter only).

Title	Action	Reason
Problem detected with the catheter	Replace Catheter	Catheter proximal sensor disconnected (vCLAS Catheter only).
Console Error Detected	If the problem persists, contact technical support	Liquid N2 latch open.

Table 13 – List of Safe Abort Pop-ups**8.4. Working with System Informational Messages and Error Pop-ups**


If an Informational Message or Error Pop-up appears, perform the following:

1. Carefully read and follow the instructions in the window.
2. If necessary, call Adagio Medical Technical Support for help.
3. If the problem is not corrected, the message will automatically reappear.

8.5. Troubleshooting the Monitor

If there is a monitor failure during a procedure, cycle power using the On/Off NUC Power Button on the rear of the monitor cover. If monitor still does not function properly, remove the rear cover on the back of the monitor. Verify if power, video, and USB cables are properly connected between the Console and Touch Screen monitor. If screen and menus are still visible, restart system by shutting down console with On/Off NUC Power Button then Main AC Power switch. If Touch Screen still fails contact Adagio Medical, Inc. If menus are no longer visible push the physical E-Stop button or On/Off NUC Power Button to power down the console. Discontinue procedure. Keep in mind catheter shouldn't be moved if an ice ball is still on the Cryoablation Element, sufficient time must elapse to allow the ice to melt.

If the monitor flickers or remains black but the console appears to be powered on, verify that the Main AC Power Switch on the top of the monitor is turned on. If necessary, cycle the button off and on to reset the touch screen. If problem persists, contact Adagio Medical, Inc.

*Figure 112 – Monitor Buttons***8.6. Console Failure during Procedure**

If the computer is still live, properly shut down the system. If not, proceed to power down the console by pressing the On/Off NUC Power Button. Discontinue the procedure. Keep in mind the catheter should not be moved or manipulated if ice is still present on the Cryoablation Element. Ensure sufficient time has passed to allow the ice to melt according to the Catheter IFU.

8.7. Sensor Failure

Recommend sensor replacement or procedure termination if possible, otherwise be extremely sensitive to patient safety and condition.

8.8. Catheter Failure

Replace Catheter. If failure persists contact Adagio Medical, Inc.

8.9. Software Failure

Restart system by shutting down the console with On/Off NUC Power Button and then turn off the Main AC Power Switch. If Software still fails, contact Adagio Medical, Inc.

9. CLEANING



Use 70% isopropyl alcohol to wipe the exterior surfaces of the console and other non-sterile components before and after use. Use a screen cleaner for the touch screen.

DO NOT clean or reuse sterile components.

Note – Discard all used catheters and sterile components in accordance with hospital procedures and the catheters' Instructions for Use.


10. SYMBOLS

The following symbols appear in the Console documentation and labels. These internationally recognized symbols are defined by the International Electrotechnical Commission, IEC 60417 and ISO 7000, and ISO 15223-1.

Symbols	Description of Symbol	Applicable Standard and Reference
	Indicates the medical device manufacturer, as defined	ISO 15223-1 5.1.1 ISO 7000 - 3082
	Indicates the authorized representative in the European Community	ISO 15223-1 5.1.2
	Indicates the date when the medical device was manufactured	ISO 15223-1 5.1.3 ISO 7000 - 2497
	Indicates the manufacturer's batch code so that the batch or lot can be identified	ISO 15223-1 5.1.5 ISO 7000 - 2492
	Indicates the manufacturer's catalogue number so that the medical device can be identified	ISO 15223-1 5.1.6 ISO 7000 - 2493
	Indicates the manufacturer's serial number so that a specific medical device can be identified	ISO 15223-1 5.1.7 ISO 7000 - 2498
	Indicates a medical device that should not be used if the package has been damaged or opened	ISO 15223-1 5.2.8 ISO 7000 - 2606
	Indicates a medical device that needs to be protected from moisture	ISO 15223-1 5.3.4 ISO 7000 - 0626
	Indicates the temperature limits to which the medical device can be safely exposed	ISO 15223-1 5.3.7 ISO 7000 - 0632
	Indicates the range of humidity to which the medical device can be safely exposed	ISO 15223-1 5.3.8 ISO 7000 - 2620
	Indicates the range of atmospheric pressure to which the medical device can be safely exposed	ISO 15223-1 5.3.9 ISO 7000 - 2621
	Indicates the need for the user to consult the instructions for use	ISO 15223-1 5.4.3 ISO 7000 - 1641
	General warning sign	ISO 7010-W001
	Non-ionizing Electromagnetic Radiation	IEC 60147/ISO 7000-5140
	Defibrillation Proof Type CF Applied Part	IEC 60417/ISO 7000-5336
	Pushing prohibited	ISO 7010-P017
	Voltage AC	N/A
50-60 Hz	Frequency (50-60 Hertz)	N/A

Symbols	Description of Symbol	Applicable Standard and Reference
3.6 A	Amperage (3.6 Amperes, 120 VAC)	N/A
2.0 – 1.8 A	Amperage (2.0 – 1.8 Amperes, 220 – 240 VAC)	N/A
EU 	Fuse rating	N/A
	Must read Operator's Manual	N/A
	Indicates the item is a medical device	ISO 15223-1 5.7.7
	Adagio Cryoablation Catheter is for exclusive use with Adagio Cryoablation Console	N/A
	CE Mark	N/A
	Indicates the entity importing the <i>medical device</i> into the locale	ISO 15223-1 5.1.8 ISO 7000 - 3725
	Indicates the entity distributing the <i>medical device</i> into the locale	ISO 15223-1 5.1.9 ISO 7000 - 3724

11. LABELING



8	Fill Port Label		16	Catheter Connection Port Label	
			17	Adagio Medical Company Logo	

Figure 113 – Console Front, Side, and Top (Fill Port) with Labels

1	Product Label			9	Main Power Label	
2	Voltage Setting Label	EU: 220-240VAC	US & CAN: 120VAC	10	E-Stop Label	
3	Weight Label			11	Power Button Label	
4	Caution Label: "CAUTION! DO NOT DISCONNECT WHILE PRESSURIZED"			12	HDMI Port Label	
5	Do Not Push, Read Label and Operator's Manual			13	USB Port Label	
6	Tank Pressure Transducer Label			14	PMAT Port Label	
7	High Pressure N2 Label			15	Grounding Port Label	
18	Intertek Certification Label					

Figure 114 – Console Rear with Labels

12. INCIDENT REPORTING

Serious incidents in relation to the device should be reported to the Manufacturer and the Competent Authority of the users’ country.

Manufacturer

Adagio Medical, Inc.
26051 Merit Circle, Suite 102
Laguna Hills, CA 92653,
USA
☎ 949-348-1188
📠 949-348-1866

European Authorized Representative

QualRep Services B.V.
Utrechtseweg 310 – Bldg B42
NL-6812 AR Arnhem
The Netherlands
Phone: +31 20-788-2630

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The screenshot shows the 'Patient Information' screen on the Adagio Cryoablation Console. It features a 'BACK' button at the top left. Below it are input fields for 'Case Number *', 'Physician Name', 'Center Name', and 'Study Name'. There is also a 'Study Date' field with a date picker. At the bottom, there is a 'Comments' text area and a 'CONTINUE' button. A legend indicates that fields with an asterisk are required.

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Adagio Medical GmbH