
vCLAS™ Cryoablation Catheter Instructions for Use

Model:

vCLAS™ Cryoablation Catheter

Basic UDI-DI- 0850037705109-0017-001D2

Confidentiality Statement

The information contained in this document is the property of Adagio Medical, Inc. It is therefore provided to you in confidence. It is understood that this information will not be disclosed to others without written authorization from the Adagio Medical Inc., Clinical and/or Regulatory Affairs.

1. DEVICE AND SYSTEM DESCRIPTION

The Adagio Medical vCLAS™ Cryoablation Catheter (VT Cryoablation Catheter or Catheter) is a single use device intended for percutaneous delivery and endocardial ablation of ventricular substrate in the treatment of ventricular tachycardia. It is designed to be used together with Adagio Medical Cryoablation Console (Console and software) and related components as a part of Adagio's VT Cryoablation System.

The Catheter connects to the Cryoablation Console via Proximal Connector (**Figure 1**) and Console's Connection Port, creating a closed loop that circulates cryogenic near-critical nitrogen (CN2) from the Console to the 15 mm long Cryoablation Element of the Catheter. This enables the Cryoablation Element to reach cryogenic temperatures and to create therapeutic, transmural linear and focal lesions. The Cryoablation Element temperature is continuously monitored by two internal thermocouples and displayed on the Cryoablation Console.

The Catheter is compatible with approved 10 F (or greater) steerable and fixed-curve guide sheaths (see **Section 9.1**). Vascular access is obtained by standard femoral vein and/or femoral artery puncture.

The Cryoablation (Freezing) Element of the Catheter contains eight (8) 1 mm electrodes distributed evenly with 1 mm spacing (**Figure 1**). The electrodes provide EGM recording and pacing capability when no ice formation is present. The EGM monitoring and pacing are accomplished via standalone Electrophysiology (EP) Recording System. The EP Recording System is a standalone medical device used for recording and monitoring of electrocardiograms (ECG) to map and diagnose ventricular tachycardia. The Catheter connects to the EP recording system via Distal Electrode Connector and standard adapter cable. The EP Recording systems use universal junction box and adapter cables. The Catheter is compatible with any junction box of commercially available EP Recording Systems.

The distal portion of the Catheter (**Figure 1**, Articulation Segment) supports bi-directional deflection of Cryoablation Element by 180°, as controlled by the actuator and turning knob within Catheter's Handle Assembly. The turning knob adjusts the resistance of the deflection up to positional lock at any desired deflection angle.


Figure 1: vCLAS™ Cryoablation Catheter

1.1. Contents of Package

The Catheter is supplied sterile. The contents include:

- One (1) vCLAS™ Cryoablation Catheter
- Product information/documentation

1.2. User Information

- Ablation procedures using the vCLAS™ Cryoablation Catheter should only be performed in the fully equipped electrophysiology laboratories, by trained healthcare professionals specialized in cardiac electrophysiology and trained in catheter manipulations and cardiac access via transeptal and/or retrograde access.
- The Catheter shall only be used by healthcare professionals who have been trained by Adagio Medical personnel regarding the proper use and advised of any safety recommendations of this Catheter.

1.3. Patient Population

The Adagio Medical Inc. VT Cryoablation System (Catheter and Console) is used in patients suffering from recurrent monomorphic ventricular tachycardia. Reference section Contraindications for more details.


1.4. Clinical Benefits

The intended clinical benefit of the Adagio Medical VT Cryoablation System is to prevent the reoccurrence of sustained VT episodes or appropriate ICD Intervention. In the CryoCureVT Clinical Study (NCT 04893317), acute success was reported at 94% non-inducible VT at the end of the ablation procedure. At 6-month follow up, 63% of the study population remained free from VT recurrence or appropriate ICD intervention after a single cryoablation procedure.

A Summary of Safety and Clinical Performance of the device is available in Eudamed (<https://ec.europa.eu/tools/eudamed/#/screen/home>).

2. LABELING AND/OR PACKAGING SYMBOLS

Symbols	Description of Symbol	Applicable Standard and Reference
	Indicates the medical device manufacturer, as defined	ISO 15223-1 5.1.1 ISO 7000 - 3082
	Indicates the authorized representative in the European Community	ISO 15223-1 5.1.2
	Indicates the entity importing the <i>medical device</i> into the locale	ISO 15223-1 5.1.8 ISO 7000 - 3725
	Indicates the entity distributing the <i>medical device</i> into the locale	ISO 15223-1 5.1.9 ISO 7000 - 3724
	Indicates the date after which the medical device is not to be used	ISO 15223-1 5.1.4 ISO 7000 - 2607
	Indicates the manufacturer's batch code so that the batch or lot can be identified	ISO 15223-1 5.1.5 ISO 7000 - 2492
	Indicates the manufacturer's catalogue number so that the medical device can be identified	ISO 15223-1 5.1.6 ISO 7000 - 2493
	Indicates the manufacturer's serial number so that a specific medical device can be identified	ISO 15223-1 5.1.7 ISO 7000 - 2498
	Indicates the medical device (vCLAS™ Catheter) that has been sterilized using ethylene oxide	ISO 15223-1 5.2.3 ISO 7000 - 2501
	Indicates the medical device that should not be used if the package has been damaged or opened	ISO 15223-1 5.2.8 ISO 7000 - 2606
	Indicates the medical device that needs protection from light sources	ISO 15223-1 5.3.2 ISO 7000 - 0624
	Indicates the medical device that needs to be protected from moisture	ISO 15223-1 5.3.4 ISO 7000 - 0626
	Indicates the temperature limits to which the medical device can be safely exposed	ISO 15223-1 5.3.7 ISO 7000 - 0632
	Indicates the medical device that is intended for one use, or for use on a single patient during a single procedure	ISO 15223-1 5.4.2 ISO 7000 - 1051
	Indicates the medical device that is not to be re-sterilized. (The vCLAS™ Catheter should not be re-sterilized.)	ISO 15223-1 5.2.6 ISO 7000 - 2608
	Indicates the need for the user to consult the instructions for use	ISO 15223-1 5.4.3 ISO 7000 - 1641

Symbols	Description of Symbol	Applicable Standard and Reference
	Indicates the need for the user to consult the instructions for use for important cautionary information such as warnings and precautions that cannot, for a variety of reasons, be presented on the medical device itself.	ISO 15223-1 5.4.4 ISO 7000 – 0434A
	Indicates the medical device that is non-pyrogenic	ISO 15223-1 5.6.3 ISO 7000 - 2427
	vCLAS™ Cryoablation Catheter is for exclusive use with Adagio Cryoablation Console	N/A
	Indicates the item is a medical device	ISO 15223-1 5.7.7
	Indicates the original medical device information has undergone a translation which supplements or replaces the original information	ISO 15223-1 5.7.8 2019-11-01
	Indicates the medical device that contains substance (Cobalt) that can be carcinogenic and reprotoxic	ISO 15223-1 5.4.10 ISO 7000 - 3723
	CE mark	N/A

3. INTENDED USE

The vCLAS™ Cryoablation Catheter, a component of the Adagio Medical Inc. VT Cryoablation System, is intended to create transmural linear and focal lesions to cardiac structures within the endocardium of the ventricles with cryoablation.

4. INDICATIONS FOR USE

The Adagio Medical Inc. VT Cryoablation System (Catheter and Console) is indicated for the treatment of monomorphic ventricular tachycardia by ablation of arrhythmogenic tissue that drives and maintains these arrhythmias.

5. CONTRAINDICATIONS

The Adagio Medical Inc. VT Cryoablation System (Catheter and Console) is contraindicated as follows:

- In patients with a systemic infection
- In patients with intracardiac thrombus
- In patients with cryoglobulinemia
- In pediatric patients (below the age of 18)
- In pregnant or breastfeeding patients
- In patients with hypercoagulopathy or an inability to tolerate anticoagulation therapy during an electrophysiology procedure
- In patients with an interatrial baffle or patch as the transseptal puncture could fail to close
- In patients with a contraindication to an invasive electrophysiology procedure where insertion or manipulation of a catheter in the cardiac chambers is deemed unsafe

6. WARNINGS AND PRECAUTIONS

- **Anticoagulation** – Administer appropriate peri-procedural anticoagulation therapy per standard of care for patients undergoing cardiac ablation therapy in left and right ventricles. Administer anticoagulation therapy during and post-procedure according to local institution standards to minimize bleeding and thrombotic complications.
- **Approved Guide Sheaths** – Use only approved guide sheaths. Refer to **Section 9.1** for the brands and models of the approved sheaths to be used with the Catheter.
 - The use of an unapproved guide sheath may cause damage to the Catheter and to the sheath and may lead to serious patient injury.
 - Follow the guide sheath manufacturer's IFU to deploy the sheath. For the Mullins sheath, always use the supplied dilator in the sheath lumen before entering venous/arterial access to avoid sheath damage.

- **Biohazard Disposal** – Dispose of the Catheter in accordance with hospital procedures and local regulatory and biohazard standards.
- **Catheter Manipulation** – Extreme care must be utilized during catheter manipulation to avoid Catheter damage, cardiac damage and/or perforation. The Catheter should be deployed in the middle of the cardiac chamber.

Do not use excessive force to advance, withdraw, deflect, or rotate the Catheter, especially if resistance is encountered. Excessive force may lead to catheter damage, including kink or cryogen leakage.

- **Carcinogenic, Mutagenic, Reprotoxic (CMR) Substance** – The device is known to contain Cobalt (CAS #7440-48-4), a CMR substance. All materials with direct patient contact were tested and met Biocompatibility requirements.
- **Cryo-adhesion** – During cryoablation, the cryoablation element adheres to the cardiac tissue and may adhere to the guide sheath. Do not attempt to manipulate the Catheter during or after the cryoablation cycle until sufficient evidence for lack of adhesion to tissue or guide sheath is observed by ensuring the Catheter Distal Temperature (T_D) reaches a minimum of 30°C, as indicated by color change from blue to white, the Thaw cycle as defined by the console software is complete, and interference seen on the EGMs during the freeze disappears. Do not rotate, pull, remove, or otherwise manipulate a Catheter that resists movement.
- **Device Integrity** – Do not use Catheter if visual inspection shows bumps, kinks, or other damage. If there is a suspicion the Catheter may be damaged or kinked while in patient, it should be carefully removed, visually inspected, and a FULL FUNCTIONAL TEST (Console Operator's Manual **Section 4.4 and 4.5.2**) should be performed outside of the patient. If any damage is confirmed, the Catheter shall not be reinserted in the patient and must be returned to Adagio Medical, Inc.
- **Device Sterility** – Catheter is supplied sterile. Inspect sterile packaging (all seals) for any signs of damage such as holes or compromised seal integrity prior to opening and removal of the device. If damaged, do not use it and return to Adagio Medical, Inc.
 - The catheter is for single use only. Do not reuse, reprocess, or re-sterilize, as those will compromise device integrity and function and may lead to patient illness, injury, or death.
- **Embolism Risk** – Introducing catheters and sheaths into the circulatory system may result in air emboli. Use proper endovascular techniques to minimize the risk of embolization, including flushing the sheath prior to catheter insertion.
- **Interaction With Implantable Cardiac Defibrillators (ICDs) and ICD Leads** – VT patients are often implanted with ICDs as protections from Sudden Cardiac Death. Disable ICD operation prior to performing cryoablation and use temporary devices to support pacing and defibrillation during the procedure. Do not perform cryoablation when Cryoablation Element is in direct contact with ICD lead(s).
- **Interaction With Other Devices** – Use caution when there are diagnostic catheter(s) in the cardiac chamber in addition to the VT Catheter. No freeze should be initiated unless it is confirmed that other catheters are not in contact with the Catheter's cryoablation element. To avoid entangling the catheters, it is recommended to manipulate different catheters sequentially, keeping them as far apart as possible.
- **Package Integrity** – VT Catheter is supplied individually packaged within a pouch and a box. Package contents are listed on the pouch and box labels. Do not use it if labeling is illegible or incomplete. Check expiration date shown on the label. Do not use the Catheter beyond the expiration date.
- **Pressurized Cryogen** – Catheter contains pressurized NCN cryogen during operation.
 - The Catheter reaches cryogenic temperatures. Do not touch the cryoablation element when cryogenic energy is delivered by the Cryoablation Console as operator injury may occur.
 - The application of cryogenic temperature results in permanent damage to cardiac tissues. Ensure the cryoablation element is in contact only with the target tissue prior to initiating energy delivery.
 - If the Catheter develops an external leak, patient injury (such as air/gas embolism) or death may occur. **CONTINUOUSLY** monitor the catheter proximal temperature on the Cryoablation Console display during cryogen delivery. If the catheter's proximal temperature (T_P) decreases to below -10°C (indicating potential cryogen leak), **IMMEDIATELY DISCONNECT and REMOVE** the Catheter.
- **Real-time Cardiac Imaging** – All Catheter maneuvers within the patient (advancement, retraction, deflection, manipulations) must be performed under fluoroscopic guidance or using other real-time imaging modalities capable to accurately identify anatomic and relative position of the devices with sufficient accuracy. Compatible imaging modalities are specified in **Section 9.3**.
- **Repeat Insertion of the Catheter** – Repeat insertions of the catheter into the guide sheath increases the risk of catheter damage. Visually inspect the catheter and perform a FULL FUNCTIONAL TEST before every reinsertion.
- **Transitory Cardiac Ischemic Events** – Cryoablation carries small risk of transitory cardiac ischemic events. Monitor for signs of transitory cardiac ischemia using 12-lead ECG and/or internal cardiac electrograms. Stop energy delivery if clinically meaningful signs of ischemia are detected and allow to fully thaw. Reassess ablation location using imaging modalities and Catheter electrograms displayed on EP System and relocate the Catheter as necessary.

- **X-rays Exposure** – The use of fluoroscopy during electrophysiology procedures may present the potential for significant x-ray exposure. Attention should be given to potential radiation exposure associated with the cryoablation procedure, taking steps to minimize exposure.

7. POTENTIAL ADVERSE EVENTS AND RISKS

The use of the Adagio Medical VT Cryoablation System does not expose the patient to significant additional risks beyond the currently available ablation technologies including RF catheters. Risks and potential adverse events of the VT Cryoablation System include but are not limited to:

- Acute respiratory distress syndrome (ARDS)
- Acute periprocedural hemodynamic/cardiogenic shock
- Air/gas embolism
- Allergic reaction/anaphylaxis
- Anemia
- Valve damage such as Aortic and/or Mitral, Tricuspid
- Aorto-right atrial fistula
- Anesthesia reaction
- Arrhythmia recurrence
- Arteriovenous (AV) fistula
- Heart block
- Bundle branch block
- Cardiac arrest
- Cardiac perforation/tamponade
- Cardiac thromboembolism
- Cerebral infarct (hemorrhagic or thromboembolic)
- Chest pain/discomfort
- Congestive heart failure/ pulmonary edema
- Coronary artery damage
- Death
- Dislodgement of pacemaker
- Electric shock
- Elevated cardiac enzymes
- Endocarditis
- Femoral nerve injury/palsy
- Femoral artery occlusion
- Femoral vein occlusion
- Fever
- Fluid overload
- Groin hematoma/ecchymosis/vascular complications
- Hypotension/Hypertension
- Hypothermia
- Infection/sepsis
- Left ventricular pseudoaneurysm
- Major bleeding, requiring surgery or transfusion
- Myocardial infarction
- Obstruction or perforation or damage to vascular system
- Pericardial effusion
- Pericarditis
- Pleural effusion
- Pneumonia
- Prolonged procedure
- Prolonged recovery
- Pulmonary bronchi thermal injury
- Radiation injury
- Respiratory depression
- Retroperitoneal bleed
- Skin burns
- ST segment elevation
- Thromboembolism
- Tissue damage related to cryo injury
- Transient ischemic attack (TIA)
- Vasovagal reactions

8. INSTRUCTIONS FOR USE

The following handling instructions shall be performed after the Cryoablation Console has been unpacked and prepared for use according to the console manual (Adagio P/N: 108-0063-001).

8.1. Preparing the vCLAS™ Cryoablation Catheter

- 8.1.1. Inspect the Catheter package for labeling and confirm the catheter type. Check the integrity of the labeling and package sterility. Do not use if compromised.
- 8.1.2. Using sterile technique, remove the catheter from its package and place in a sterile area.
- 8.1.3. Visually inspect the catheter for damage and overall condition. If impaired, do not use it. Return to Adagio Medical, Inc.

8.2. Connect Catheter to Console & Test

- 8.2.1. Prepare the VT Cryoablation Console per Operator's Manual (P/N 108-0063-001) **Section 4.3**.
- 8.2.2. Using Console GUI initiate Catheter connection per **Section 4.4** of the Console Operator's Manual (P/N 108-0063-001). Maintaining the sterility of the catheter, align and insert the Catheter Proximal Connector into Console's Connection Port. An audible click will indicate that the Console has secured the Catheter Connector in the Port.
- 8.2.3. The Console will test the Catheter's vacuum and perform self-tests per **Steps 7 to 9** of the Operator's Manual (P/N 108-0063-001) **Section 4.4**. During self-tests the Catheter Cryoablation Element may briefly reach cryogenic temperatures. Ensure that Cryoablation Element does not come into contact with the body or any other surfaces during self-tests.

- 8.2.4. **FULL FUNCTIONAL TEST** - Immerse the Cryoablation Element of the Catheter into a container with the sterile saline solution and initiate **FULL FUNCTIONAL TEST** per **Step 10** of Operator's Manual (P/N 108-0063-001) **Section 4.4**. Observe ice forming on the Cryoablation Element.
- Verify that **NO gas bubbles** are released from the cryoablation element, consistent ice formation, and the Catheter Proximal Temperature (T_p) on Console GUI screen remains above -10°C .
 - The lack of ice formation, **observed bubbles**, or Catheter Proximal Temperature (T_p) below -10°C represent the failure of FULL FUNCTIONAL TEST. The catheter **MUST** be replaced.
 - Upon completion of the FULL FUNCTIONAL TEST allow ice to melt completely off the Cryoablation element. If the Cryoablation Element adheres to the container wall during freeze do not apply force to free it up – allow sufficient melting of the ice to release the Catheter prior to further manipulation.
- 8.2.5. If the tests described in sections above are successful, the Catheter is ready for use. Connect the Catheter's Electrical Connector to the male connector on the compatible EP System's Adapter Cable (see section 9.2 for the list of the compatible cables). Connect the other end of the Adapter Cable into the EP System's Catheter Interface; note that the Freezing Tip contains all of the electrodes on Adagio's Catheter and is connected to pin # 1 – pin # 8 on the junction box.
- 8.2.6. If the tests are NOT successful, disconnect the Catheter following the steps in **Section 4.6** of the Operator's Manual and return the Catheter to Adagio Medical, Inc. Prepare a new Catheter beginning from **Section 8.1** above and repeat test sequence.

8.3. Target Chamber Access And Diagnostic Maneuvers

- 8.3.1. Using an aseptic technique, create an appropriate vascular access in the femoral veins and/or femoral arteries sufficient for direct (for RV), transseptal or retrograde (for LV) positioning of Cryoablation Catheter and diagnostic catheters in the cardiac chambers of interest. Initiate anticoagulation protocol per institutional standard of care. Using fluoroscopic and/or intracardiac echo (ICE) guidance perform transseptal puncture, as necessary.
- 8.3.2. Place appropriate diagnostic catheters for the ablation procedure.
- 8.3.3. Using fluoroscopic guidance place Adagio approved ≥ 10 F guide sheaths (reference **Section 9.1** for specific models), insert and advance the Catheter to the cardiac chamber of interest. Minimize deflection of the transeptal guide sheath prior to advancing the Catheter and ensure the Catheter exits the sheath in the middle of the cardiac chamber. Use Catheter deflection mechanism to facilitate proper positioning and crossing of cardiac valves, per institutional standard of care.
- 8.3.4. Clip the proximal portion of the catheter (gray proximal shaft) to the table or the drape to prevent accidental or unintentional movement.
- 8.3.5. Perform required diagnostic maneuvers, inclusive of VT induction, entrainment, pace, activation and/or substrate mapping. If using Electroanatomic Mapping System (EAMS), ensure that the Catheter connections to the EP Station junction box are routed to the EAMS inputs and that the Cryoablation Element of the Catheter is properly rendered within the System.
- To check the EGM signal with the Catheter, advance it to the cardiac chamber of interest and manipulate using standard endovascular techniques and catheter/sheath deflections, as appropriate, guided by fluoroscopy and/or other imaging modalities.

8.4. Delivery of Cryoablation

- 8.4.1. To perform cardiac ablation, advance the Catheter to the target tissue, using standard endovascular techniques and catheter/sheath deflections, as appropriate, guided by fluoroscopy and/or other imaging modalities. Do not rotate, push/pull or otherwise manipulate a Catheter that resists movement.
- For linear lesions: to maximize the delivery of cryogenic energy, position the Catheter Cryoablation Element in parallel to the target tissue ensuring maximum tissue contact along the entire length of Cryoablation Element.
 - For focal lesions: ensure that the most distal part of Catheter Cryoablation Element (at least 25% of the total length) is in contact with the target tissues.
- 8.4.2. Ensure the Console is in the Procedure Ready State per Operator's Manual **Section 4.4** and initiate cryoablation cycle. Successful cryoenergy applications will be administered in cycles. Based on tissue thickness, one complete cycle includes an initial freeze ranging from 30 seconds to 5 minutes. This is followed by a thaw of one minute (minimum) and a bonus freeze which shares the same duration as the initial freeze to complete one full cycle (freeze-thaw-freeze cycle).
- Prior to initiation of the freeze, ensure that no other catheter is in direct contact with Cryoablation Element.
 - Use cryomapping if there is a clinical suspicion of catheter proximity to the phrenic nerve or cryoablation-sensitive parts of the cardiac electrical conduction system.
 - During cryoablation, the cryoablation element adheres to the cardiac tissue and may adhere to the guide sheath. Do not attempt to manipulate the Catheter during or after cryoablation cycle until sufficient evidence for lack of adhesion to tissue or guide sheath is observed by ensuring the Catheter Distal Temperature (T_D) reaches a minimum of 30°C , as indicated by color change from blue to white, the Thaw cycle as defined by the console software is complete, and interference seen on the EGMs during the freeze disappears. Do not rotate, pull, deflect, or otherwise manipulate a Catheter that resists movement.

- 8.4.3. During cryoablation monitor for signs of transitory cardiac ischemia using 12-lead ECG and/or internal cardiac electrograms. Stop energy delivery if clinically meaningful signs of ischemia are detected and allow to fully thaw. Reassess ablation location using imaging modalities and Catheter electrograms displayed on EP System and relocate the Catheter as necessary.
- 8.4.4. Repeat **Steps 8.4.1. to 8.4.3.** above until all target tissue is ablated and all diagnostic maneuvers are completed.

8.5. Catheter Removal and (Optional) Reinsertion

- 8.5.1. Ensure the Catheter is fully thawed ensuring the Catheter Distal Temperature (T_D) reaches a minimum of 30°C, as indicated by color change from blue to white, the Thaw cycle as defined by the console software is complete, and interference seen on the EGMs during the freeze disappears. Retract Cryoablation Element under fluoroscopic guidance into the guide sheath while the Catheter and guide sheath are not deflected to ensure low retraction force. Do not rotate, pull, deflect, or otherwise manipulate a Catheter that resists movement.
- 8.5.2. Carefully withdraw the Catheter from the guide sheath and place it in the sterile area.
- 8.5.3. If reinsertion of the Catheter is required, visually inspect the Catheter per **Section 8.1** and perform FULL FUNCTIONAL TEST per **Steps 8.2.4.** above prior to inserting back into the guide sheath.
- 8.5.4. To end the procedure and fully disconnect the Catheter follow the Console Operator's Manual **Section 4.6.**

9. DEVICE COMPATIBILITY

9.1. Adagio Approved Guide Sheaths

The use of a guide sheath unapproved by Adagio may damage the catheter and the sheath. Use only approved guide sheaths:

- Oscor Destino Reach (10 F)
- Abbott Fast-Cath Mullins (10 F)

9.2. Compatible Adapter Cables

- Boston Scientific 5454S
- Biosense Webster C6MRMST10SA

9.3. Imaging Modalities

A variety of imaging modalities can be used with the Adagio Medical Inc. VT Cryoablation System (Catheter and Console). These modalities are chosen by the physician based on preference. The following are used but are not limited to:

- Fluoroscopy
- Intracardiac echocardiography (ICE)
- Electroanatomical Mapping (EAM) System (e.g., Abbott Ensite™ NavX™, Boston Scientific Rhythmia HDx, and Biosense Webster Carto®)

10. SPECIFICATIONS

Catheter Shaft Diameter	9F
Catheter Usable Length	115 cm
Freezing Element / Tip Length	15 mm
Compatible Sheath	≥ 10 F
Electrode Size	1 mm
Electrodes	Freeze Element: 8
Electrode Configuration	Freeze Element 1-1-1 repeated
Cryoablation Performance	≤ -100°C achieved within 30 seconds
Environmental Parameters	Non-pyrogenic Single Use Sterilization: EO

11. STORAGE

Store the vCLAS™ Catheter in a cool, dry location. Avoid excessive heat and cold during transport. Recommended temperature for storage is +15°C to +22°C. Recommended temperature for transportation is -5°C to +55°C.

12. DISPOSAL

Used products are contaminated and must be handled and disposed of as contaminated hospital waste.

13. INCIDENT REPORTING

Serious incidents in relation to the device should be reported to the Manufacturer and the Competent Authority of the users' country.

Manufacturer

Adagio Medical, Inc.
26051 Merit Circle, Suite 102
Laguna Hills, CA 92653,
USA

 949-348-1188

 949-348-1866

European Authorized Representative

QualRep Services B.V.
Utrechtseweg 310 – Bldg B42
NL-6812 AR Arnhem
The Netherlands
Phone: +31 20-788-2630

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Adagio Medical GmbH